


| Re-Order # | Size | Unit |
|------------|-----------|------|
| 911585 | 24" x 72" | Each |
| 911590 | 30" x 72" | Each |



Beveled Low Profile Edges

Purpose

Skil-Care EZ Landing Fall Mats are designed to help reduce the impact from bed related falls. The one-piece electrometric foam surface is beveled to provide a low-profile edge to reduce tripping hazard. The anti-microbial agent protects the mat from mold, mildew and odor causing bacteria. The low-profile surface allows resident to stand, walk or wheel over the mat without difficulty. The mat comes with a non-skid top and a slip resistant bottom.

Easy to Use

This mat is available in two widths and the choice of width should be made by a trained professional.

1. Place the mat in front of the bed spaced evenly between both legs and parallel with the bed (mattress).
2. Assure that the mat is placed with the beveled side up.
3. Although this product is durable, rolling the bed wheels over the mat should be avoided.
4. The mat could be stored flat or rolled loosely.

Easy Care

The EZ Landing II fall mat could be wiped clean with a spray and wipe cleanser. If disinfected, avoid using chlorine (bleach) chemicals.

Turn product insert over for more of Skil-Care's Fall Mat products.

Alternative **Skil-Care**™ CORPORATION **Fall Mats.....**

Beveled Edges


Roll-On Bedside Fall Mat

This mat protects patients who are high risk of falling from the bed and incorporates trip resistant beveled edges.

- High density shock absorbing packaging foam
- Step on stability
- Slip resistant bottom
- Heat-sealed vinyl
- Easy storage strap for folding or securing to the footboard



| Reorder# | Description | Size | Unit |
|----------|-----------------|------------------|------|
| 911548 | Roll-On Bedside | 1.5" x 26" x 68" | Each |
| 911549 | Roll-On Bedside | 1.5" x 36" x 68" | Each |



Beveled Edge



SafeSide Fall Mat

A fall mat that resists movement on the floor and provides trip resistant beveled edges.

- High density shock absorbing bonded foam
- Step on stability
- Slip resistant bottom
- Heat-sealed vinyl



| Reorder# | Description | Size | Unit |
|----------|----------------|----------------|------|
| 911572 | Bi-Folding Mat | 1" x 24" x 68" | Each |
| 911577 | Bi-Folding Mat | 1" x 36" x 68" | Each |

Non-Beveled Edges

Raised Height Soft Fall Bedside Mat

For maximum fall protection this is the optimum fall mat.

- High density polyethylene foam
- Slip resistant bottom
- Heat-sealed vinyl
- Easy storage strap for folding or securing to the footboard



| Reorder# | Description | Size | Unit |
|----------|-------------|----------------|------|
| 911546 | Folding Mat | 2" x 24" x 68" | Each |
| 911536 | Folding Mat | 4" x 24" x 68" | Each |
| 911545 | Folding Mat | 2" x 36" x 68" | Each |
| 911535 | Folding Mat | 4" x 36" x 68" | Each |

| Reorder# | Description | Size | Unit |
|----------|-----------------|----------------|------|
| 911540 | Non-Folding Mat | 2" x 36" x 68" | Each |
| 911530 | Non-Folding Mat | 4" x 36" x 68" | Each |



Tri-Fold Bedside Fall Mat

Designed for easy storage

- High density shock absorbing packaging foam
- Step on stability
- Slip resistant bottom
- Heat-sealed vinyl
- Velcro straps for easy folding and storage



| Reorder# | Description | Size | Unit |
|----------|-----------------|----------------|------|
| 911547 | Tri-Folding Mat | 1" x 26" x 68" | Each |



Visit us online at www.skil-care.com for our full line of nursing and therapy products.