

Case Cart

Model: 9203-C

Medical grade case carts constructed of aesthetically pleasing melamine material to create a warm medical environment. Designed to meet all of your practical needs with their various compartments and drawers



Features:


- Constructed of medical grade Melamine material
- Wide cabinet with one shelf
- Drawer can be on bottom or top (top is standard unless specified differently)
- Extender Shelf: 14"
- 5" locking casters
- 1 Year Limited Warranty

Overall Dimensions

34¼"H x 32½"W x 23"D



Finish Options:



Prestige Maple

Candlelight

Classic Cherry

Biltmore Cherry

Truffle