

Diagnostix™ ADSTATION™


Modular Diagnostic System

WARNING:

This unit is ONLY compatible with an Adstation transformer whose serial number is 100000 or higher. Extension module will be damaged when used with other devices.

CAUTION:

When plugging in the connecting cable, make sure the four contacts of the socket align with the four contacts of the pin strip.



Warranty:

This Diagnostix Adstation transformer and extension modules are warranted for five years from date of purchase. The warranty does not apply to damage caused by improper handling, accidents, improper use, or alterations made to the instrument by third parties. Damage caused by connection of incompatible extensions modules to Adstation transformers is not covered; please confirm compatibility of all components before installing. The warranty is only valid after the product is registered online at:

www.adctoday.com/support/warranty-registration.



American Diagnostic Corporation
55 Commerce Drive, Hauppauge, NY 11788

tel: 631-273-9600
toll free: 1-800-232-2670
fax: 631-273-9659

www.adctoday.com
email: info@adctoday.com

IB p/n 93-5699-00 rev 3 (11/29/22)

To register your product, visit us at

www.adctoday.com
and follow the links

For questions, comments, or suggestions,
call in the U.S.A. Toll Free:

1-800-ADC-2670

Printed in U.S.A.