



Sterile Antimicrobial Bandage Roll 4" x 4.1 yd

AB100

1/PK, 100 PK/CS

Product Features

Contains PHMB and benzalkonium chloride creating a broad spectrum antimicrobial effect within the dressing for gram-positive and gram-negative bacteria and fungi

May be used as a primary or secondary dressing

It has been shown to have 4 log bacterial reduction against staphylococcus aureus, escherlchia coli, candida albicans, pseudomonas aeruginosa and bacillus subtilis, streptococcus pyogenes, serratia marcescens, aspergillus niger

As a dressing the fluff roll provides quick wicking action, exceptional absorption, and excellent aeration

Attributes

Brand Name	Dukal
Size	4" x 4.1 yd
Sterile	Yes
Weight	10.58 lbs
Everyday Care	No
Premium Care	No

Packaging

Packaging	1/PK, 100 PK/CS
Item Packaging	1/PK
Case Packaging	100/CS
Total Case Quantity	100
Item Packaging Barcode	(01)10665973028092
Case Packaging Barcode	(01)20665973028099

