



accutest® CLIA Waived DRUG TEST CUPS



Just add urine and read results... it's that easy



- Accurate
- Easy-to-use
- Self-contained
- Rapid results in 5 minutes
- 24 months shelf life
- Stores at room temperature
- Built-in test controls
- Leak-proof cup may be sent to lab for confirmation

Drug Test Cup	Tests for	Catalog No.
5-Panel	AMP/COC/OPI/PCP/THC	DS86S625
5-Panel	AMP/COC/BZO/OPI(300)/THC	DS93S625
6-Panel	AMP/COC/MET/OPI/PCP/THC	DS87S625
11-Panel	AMP/MET/COC/THC/OPI/BZO/PCP/MTD/OXY/TCA/BAR	DS811
12-Panel	AMP/MET/COC/THC/OPI/BZO/PCP/MTD/OXY/TCA/BAR/MDMA	DS822
13-Panel	AMP/BAR/BZO/BUP/COC/MDMA/MET/MTD/OPI300/OXY/PCP/TCA/THC	DS823


*Customized drug screen configurations available with minimum order.

Tests for:

- Amphetamine (AMP)
- Barbiturates (BAR)
- Benzodiazepines (BZO)
- Buprenorphine (BUP)
- Cocaine (COC)
- Ecstasy (MDMA)
- Marijuana (THC)
- Methadone (MTD)
- Methamphetamine (MET)
- Opiates (OPI)
- Opiates (OPI)
- Oxycodone (OXY)
- Phencyclidine (PCP)
- Tricyclic Antidepressants (TCA)

Cutoff Levels:

- 1000 ng/ml
- 300 ng/ml
- 300 ng/ml
- 10 ng/ml
- 300 ng/ml
- 500 ng/ml
- 50 ng/ml
- 300 ng/ml
- 1000 ng/ml
- 2000 ng/ml
- 300 ng/ml
- 100 ng/ml
- 25 ng/ml
- 1000 ng/ml



JANT PHARMACAL CORPORATION
16530 Ventura Blvd., #512, Encino, CA 91436
800.676.5565 818.986.8530 Fax 818.986.0235
www.accutest.net info@accutest.net

Accutest is a registered trademark of Jant Pharmacal Corporation

