



# AliMed®

Helping people...help people™...for 38 years

## Familiar seatbelt reminds patients to stay seated

### ALIMED® SEATBELT ALARM

Helps keep patients safely seated in wheelchairs. It sounds an alarm the moment a resident releases the familiar automobile-style buckle. Alarm sounds at one of two volume levels and flashes a red light, making it easy for staff to quickly locate patient at risk. With the optional Nurse Call Connector with Y-Adapter, alarm may be set to sound at the nurses' station without interfering with the patient's nurse call button.

The durable nylon web belt easily attaches to wheelchair's back uprights. Multiple mounting grommets ensure a correct fit.

#70062 Seatbelt Alarm w/TR2

### ACCESSORIES

(SOLD SEPARATELY)

- #77788 Replacement Buckled Seatbelt (no alarm)
- #79205 12" Bariatric Seatbelt Extender
- #79990 Replacement TR2 Alarm
- #2035 Nurse Call Connector with 8' Cord and Y-Adapter
- Allows alarm to sound at nurses' station
  - Won't interfere with nurse call button



### KEY BENEFITS

- ▶ Seatbelt acts as a reminder for patient to stay seated
- ▶ Alarm sounds the moment buckle is released
- ▶ Alarm sounds and flashes to identify patient at risk
- ▶ Alarm features tamper-resistant, low-battery indicator and nurse call jack
- ▶ Easy to install



297 High Street, Dedham, MA 02026 USA  
www.AliMed.com  
E-mail: info@alimed.com

Order Toll Free 800.225.2610  
24-Hour Fax 800.437.2966  
General Info 781.329.2900