

# Avante SLS 2500

## LED Surgical Light

Provides a shadowless, natural color temperature

The Avante SLS 2500 LED Mobile Lighting System is a compact system. The shadowless lamps provide color temperature close to natural sunlight, reflecting the actual colors of human tissue.

### FEATURES

- Dimming knob located on the side of the lighthouse for quick and easy adjustments.
- Multipath light ensures that a single damaged LED in one path will not affect operation.
- The SLS 2500 is available as a ceiling mount! Ask your representative for details.
- Universal appending system provides maximal adjusting range, suitable for any operating room. Appending arm is made of nova alloy, with light-strong structure, convenient operation and accurate positioning.
- Lamp cap is totally enclosed with a streamlined design, complying with requirements of sterilizing and laminar flow cleaning in all kinds of operating rooms.



471-2023-02-04

### SPECIFICATIONS



**Illumination Intensity:**  
60,000 lux

**CRI - Color Rendering Index:**  
85 - 95

**Diameter of Light Spot:**  
7"

**Illumination Depth:**  
23.5"

**Power:**  
20 W



**Light Options:**



Wall-Mounted



Mobile



Ceiling Mount