

DEBUG BABY JOGGER



DEBUG BABY JOGGER

debug

The Debug Baby Jogger allows the parent or care giver to take their child safely across all terrains for exercise and fun.

LIGHT WEIGHT ALUMINUM TUBULAR FRAME / EASILY FOLDS FOR TRANSPORTATION AND STORAGE / FOOT BRAKE AND LOCKING MECHANISM PROVIDES STABILITY WHEN ON A GRADE / STORAGE AREA UNDER THE SEAT / #316L STAINLESS STEEL CASTERS AND AXLES / PNEUMATIC WHEELS ROLL ACROSS SAND, SNOW, SHELLS GRASS AND GRAVEL / LOW MAINTENANCE DELRIN BUSHINGS ARE IMPERVIOUS TO CORROSION / REMOVABLE HOOD AND TRAY / FITS THROUGH A STANDARD ADA DOORWAY

SPECIFICATIONS

SEAT SIZES	ENGLISH (INCH)	METRIC (CM)
Width	24	61
Depth	40	102
Height	42	107
Wheel size front	9.5	24
Wheel size rear	12	30

WEIGHTS		
Capacity	50 lbs	23 kg
Chair weight	32 lbs	15 kg

Manufactured in the USA,
Debug has a 30 day money back guarantee.

© Copyright Deming Designs Inc., 2001–2016. All rights reserved.

ORDER ONLINE WWW.DEBUGMOBILITYPRODUCTS.COM