


BAILEY

121814A

Model 364 Child Four Place Standing Table Assembly Instructions



STEP 1. FOR SHIPPING PURPOSE ONE ANGLE IS REVERSED, TURN ANGLE AROUND TO MATCH EXISTING (4) ANGLES.

STEP 2. ARRANGE STANDING BOXES AS SHOWN, ALIGN TOP TO BOXES AND MATCH ANGLES TO PREDRILLED HOLES UNDER TABLE TOP AND FASTEN DOWN WITH HARDWARE SUPPLIED.