

Sure-Safe® Bathtub Safety Rail

- Heavy-gauge, welded steel construction
- Clamp holds firmly to tub wall without marring tub surface
- Provides two gripping surfaces to facilitate both entry and exit from tub
- Color: White
- Fits tubs with walls 2 ¾" - 6 ½" thick
- Not recommended for use with fiber glass tubs
- Tool-free installation
- Maximum Weight Capacity 250 lb - EVENLY DISTRIBUTED
- Lifetime limited warranty

2032A 1 ea

