



BOVIE[®] 1250S-V, SPECIALIST|PRO HIGH FREQUENCY ELECTROSURGICAL GENERATOR

The Bovie 1250S-V, Specialist | PRO is the affordable solution for the busy veterinarian office and specialty facility. An improvement to the number one selling unit, the Aaron 1250-V, the Bovie 1250S-V has been enhanced with a modern push-button design and industry-standard bipolar connections. Capable of most monopolar and bipolar surgical procedures and packed with dependable safety features, the Bovie 1250S-V, Specialist | PRO satisfies the veterinarian's demands with precision, safety, and reliability.

- The precision of using a “pencil” over a scalpel saves time during procedures and allows veterinarians to perform more procedures per day.
- The ability for veterinarians to quickly switch from Cutting to Coagulation by simply pushing the buttons on the pencil during a procedure keeps them from having to leave the sterile field, reducing risk and saving time.
- With reusable, autoclavable pencil, 5 electrodes, and a stainless steel grounding plate, the A1250S-V is ready to go to work immediately and can reduce operating costs.
- Veterinarians can perform multiple procedures with one ESU system, enabling fast ROI and a new revenue stream with the ability to perform more procedures.



BOVIE® 1250S-V, SPECIALIST|PRO HIGH FREQUENCY ELECTROSURGICAL GENERATOR



Product code: A1250S-V

Standard Equipment

- 1 Bovie 1250S-V Electrocauter Generator
- 1 Power Cord
- 1 ESPR2 Sterile Reusable Pencil
- 5 Non-Sterile Reusable Electrodes (2 ES01R Blades, 2 A833 Needles, 1 ES20R Ball)
- 1 A1204P Metal Plate
- 1 A1254C Grounding Cord
- 1 MC-55-237-001_r0 User's Guide*

EQUIPMENT INCLUDED



*Not Pictured

PRODUCT SPECIFICATIONS

<p>Dimensions</p> <p>Height: 15.2cm (6 in.) Width: 26cm (10.25 in.) Depth: 30.5cm (12 in.) Weight: < 6.5 kg (< 14 lbs.)</p>	<p>Output Power</p> <p>Cut: 120 watts Blend: 90 watts Coagulation: 80 watts Fulguration: 40 watts Bipolar: 30 watts</p>	<p>Output Frequency</p> <p>350 – 600 kHz</p> <p>Line Voltage</p> <p>A1250S-V: 100-240 VAC</p> <p>Line Frequency</p> <p>50 – 60 Hz</p> <p>Warranty: <u>FOUR YEARS</u></p>
---	--	--



5115 Ulmerton Road, Clearwater, FL 33760 • U.S. Phone 1-800-537-2790 • Fax 1-800-323-1640
 Int'l. Phone +1-727-384-2323 • Fax +1-727-347-9144 • sales@boviemed.com
www.BovieMed.com/animalhealth.asp