From Dirty to Clean and Everything In Between... It's BrandMax Clean!®







Tri-Clean by BrandMax **Ultrasonic Cleaners**

Tri-Clean ultrasonic cleaners by BrandMax provide superior cleaning power which results in cleaner instruments. All Tri-Clean ultrasonic units come standard with heating element, cycle time controller and stainless steel basket, there is no need to pay extra for these valuable options when BrandMax provides them for free!

Take the *Tri-Clean* Challenge ...

FEATURE

Tri-Clean Competitor

BENEFIT

All stainless steel construction Stainless steel drain valve Superior acoustic pulse generators Heating element included Stainless steel basket included Works with all popular cassette systems 1 year warranty on unit & electronics / 2 year warranty on transducers

Scratch resistant, durable and easy to clean

Easy and effective draining of soiled cleaning solution

Powerful cavitations loosen bio-burden from instruments

Expensive option included for free!

Expensive option included for free!

Supports instrument management

Comprehensive protection for your investment

Model#	Cleaning Tank Capacity Gallons/Liters	l Rango	Timer Setting	Power Sonic/ Heat	Tank Dimensions LxWxH (in.)	Basket Dimensions LxWxH (in.)	Overall Dimensions LxWxH (in.)	Minimum Counter Opening W x D (in.)	Unit Weight (lbs.)	Shipping Weight (lbs.)	Shipping Dimensions LxWxH (in.)
RECESSED MODELS – SPECIFICATIONS FOR SN 3260 AND HIGHER											
U-10LHREC	2.64 gal / 10L	68-120 °F 20-60 °C	1 - 99 min.	220W / 350W	12 x 9.5 x 6 in.	10.5 x 8.25 x 4.25 in.	18.6 x 16.2 x 11.5 in.	14 x 11.5 in.	26 lbs.	34 lbs.	21 x 19 x 15 in.
U-13LHREC	3.43 gal / 13L	68-120 °F 20-60 °C	1 - 99 min.	300W / 450W	13 x 11.8 x 6 in.	11.6 x 10.6 x 4.25 in.	19.25 x 18.10 x 13 in.	15.5 x 14 in.	33 lbs.	41 lbs.	23 x 21 x 15 in.
U-19LHREC	4.94 gal / 19L	68-120 °F 20-60 °C	1 - 99 min.	300W / 450W	13 x 11.8 x 8 in.	11.6 x 10.6 x 6.25 in.	19.25 x 18.10 x 15 in.	15.5 x 14 in.	36 lbs.	44 lbs.	22 x 22 x 17 in.
U-20LHREC	5.06 gal / 20L	68-120 °F 20-60 °C	1 - 99 min.	430W / 650W	20 x 12 x 6 in.	18.85 x 11.25 x 4.25 in.	26.75 x 18.6 x 13 in.	23 x 15 in.	44 lbs.	53 lbs.	29 x 21 x 15 in.









Order Today! Call 1-877-561-0500

Buy Online: www.duralinesystems.com 39 Western Hwy, West Nyack, NY 10994

