

# PHC115G



63.75" x 23.63" x 24.13" (H x W x D)

Single chamber warming cabinet with glass door, stainless steel construction, digital thermostat and lock

## Highlights:

Stainless steel construction inside and out for durable use and easy cleaning

Digital controls conveniently located for easy adjustment and temperature viewing

Full length glass door keeps complete interior in view for easy inventory check

## Product Features:

Digital thermostat	Electronic controls are conveniently located for easy temperature adjustment
Electronic readout	Large LED temperature display ensures your staff can easily monitor the interior temperature
120 to 150°F Temperature Range	Designed and set for hospital blankets
Double pane tempered glass door	Heat-safe glass door offers a full view of stored contents, with a stainless steel door trim to highlight the interior
Professional handle	Stainless steel handle completes the professional look
Cantilevered shelving system	Easily adjustable to accommodate virtually any size
Stainless steel interior	Interior walls are constructed from stainless steel to ensure lasting durability and easy cleaning
Fully finished stainless steel cabinet	Complete stainless steel construction allows freestanding use
Factory installed lock	Keyed lock for a secure interior

## PHC115G Specifications:

Overview	
Height of Cabinet	63.75" (162 cm)
Width	23.63" (60 cm)
Depth	24.13" (61 cm)
Depth with Handle	25.25" (64 cm)
Depth with door at 90°	46.25" (117 cm)
Capacity	11.0 cu.ft. (311 L)
Door	Glass
Cabinet	Stainless Steel
US Electrical Safety	ETL
Canadian Electrical Safety	ETL-C
Amps	5.5
Voltage/Frequency	115 V AC/60 Hz
Weight	145.0 lbs. (66 kg)
Shipping Weight	150.0 lbs. (68 kg)
Parts & Labor Warranty	1 Year
UPC	761101044811
Warming Cabinet	
Door Swing	RHD
Reversible	No
Adjustable Shelves	Yes
Shelf Type	Wire
Shelf Qty	4
Thermostat Type	Digital
Interior Height	50.0" (127 cm)
Interior Width	19.5" (50 cm)
Interior Depth	18.75" (48 cm)
Dolly	Optional
Temperature Range	120 to 150°F



