

# Hate your crutches? Limited by your knee scooter?

The NEXT GENERATION of HANDS-FREE mobility is here!

Easily manage you daily activities, play with your kids, hold your baby, get up and down stairs, use public transit, work, exercise, or go shopping!

### **■ NEW & IMPROVED**



#### SELF STABILIING FOOT -

Balance is easy, safe and intuitive. Rotates 180 degrees without tools for left or right leg use. Replaceable rubber grip treads.



#### QUICK ON / OFF STRAPS -

On or off in seconds, comfortable woven straps with comfort padding. Intuitive, infinite adjustment.



#### EASY GRIP HANDLE -

Makes learning even easier, allows short trips without strapping in.



#### ERGONOMIC KNEE PLATFORM -

Contoured for all day comfort and enhanced stability. Exogrid construction for strength and stiffness.



ADJUSTABLE THIGH SUPPORTS -

Easy tool free adjustment for both thigh size and vertical alignment

"It's not about what you CAN'T do, it's about what you CAN do!"

— Erin Earlywine









■ LIFESTYLE

# **GET YOUR LIFE BACK**

EAS

EASY TO LEARN
EASY TO USE

PROVEN

AFFORDABLE

COMFORTABLE

"iWALK2.0 lets me do life" - Tim Chavez









Lightweight, durable, comfortable and inexpensive, the world's best selling hands-free crutch just got even better.

Fits heights 4'10" - 6'6", 275 lbs

(147-198 cm, 125 kg)

#### WARNING A

For lower leg injuries only. Should not be used if balance or leg strength is compromised.









## For lower leg injuries including:

Achilles injuries, foot / ankle sprains and fractures, foot ulcers, tibia / fibula fractures Most lower leg non-weight bearing injuries

## BENEFITS:

- Eliminates the pain, fatigue and inconvenience of conventional crutches
- Improves mobility during recovery period
- Easier to learn than conventional crutches
- Boosts user independence, morale and confidence
- Fully adjustable Fits 4'10" to 6'6"
- For left or right leg use reconfigure in seconds without tools
- Tool free assembly and fitting





"As a physician, I've never seen a device that does so much for the patient."

- Dr. Gretchen Lorenz, MI

"Having now used the iWALK2.0, I would never consider using crutches again. Neither should you."

"I will never be able to describe in words how much the iWALK2.0 has done for me and my family.



Easy as 1-2-3



Easy, quick, tool free



Self fit in minutes



One small step for man One giant leap for mobility



# iwalk-free.com

The world's original hands free crutch

Dealer information here

iWALKFree, Inc. 956513 C7 Mulmur, Mansfield, ON LON1MO, Canada 877 268-1466 sales@iwalk-free.com

## **LEARN MORE HERE!**

For lower leg non-weight bearing injuries



US and International Patents: www.iwalk-free.com/patent

HANDS FREE PAIN FREE



# **KNEE CRUTCH**

GET YOUR LIFE BACK





