Respiratory roscoemedical.com

Aerosol Therapy

Roscoe Bunny Pediatric Nebulizer System

Designed to reduce anxiety during treatments, the child-friendly Roscoe Bunny Nebulizer provides effective treatments with a state-of-the-art piston compressor. Delivers 43 psi at maximum pressure and 16 psi operating pressure (liter flow range 5 - 8 LPM) and weighs only 3.4 pounds. Equipped with a built-in carrying handle and a convenient nebulizer parking post. Comes complete with compressor, nebulizer, child mask, tubing, angled mouthpiece, five extra filters and a five year warranty. The Roscoe Bunny Nebulizer Compressor is also available with the TruNeb Reusable Neb Kit option.









NEB-BUNNY



Carrying Bag and Replacement Filter

Product Specifications

- Child-friendly design helps to reduce anxiety during treatments
- Provides smooth, quiet treatments
- · Equipped with a convenient nebulizer parking post
- Comes complete with compressor, nebulizer, child mask, angled mouthpiece, tubing, five extra filters and instruction manual
- TruNeb Reusable Neb Kit available purchase with the compressor or alone
- BUNNY-TRU includes TruNeb Reusable Neb Kit
- · Five year warranty
- HCPCS Code: E0570*
- HCPCS Code: A7003* for disposable nebulizer kits
- HCPCS Code: A7005* for reusable nebulizer kits
- HCPCS Code: A7013* for foam filters

Technical Specifications

- Power: I20 VAC/60Hz
- Power Consumption: 1.5 amps
- Sound Level: less than 60dBA at 1 meter
- Maximum Compressor Pressure: 25 to 42 psi
- Typical Operating Pressure: 8 to 13 psi
- Typical Operating Flow Rate: 5 ~ 8 LPM
- · Medication Capacity: 5 ml
- Particle Size: Less than 5 microns
- Operating Temperature: 50o 104oF
- Dimensions: 8.78" (L) x 5.4" (W) x 9" (H)
- Weight: 3.4 lb.

Description	Qty
Roscoe Bunny Pediatric Nebulizer System with Disposable Neb Kit	each
Roscoe Bunny Pediatric Nebulizer System with one Disposable Neb Kit and one TruNeb Reusable Neb Kit	each
Nylon Carry Bag for Roscoe Bunny Pediatric Compressor System	each
Foam Filter, ⁵ / ₈ " diameter x ½" thick	6/pk
	Roscoe Bunny Pediatric Nebulizer System with Disposable Neb Kit Roscoe Bunny Pediatric Nebulizer System with one Disposable Neb Kit and one TruNeb Reusable Neb Kit Nylon Carry Bag for Roscoe Bunny Pediatric Compressor System

