

Advanced Wound Care

Skin health is mission critical





You are the angels

We are your wings.

Cardinal Health offers a comprehensive portfolio of advanced and traditional wound care, incontinence care, skin management and negative pressure wound therapy products designed to help you maintain your patients' skin integrity and support wound healing. We believe that skin health is mission critical, and we are here with products and education to help you care for your patients and keep them comfortable — when they need it most.

Start with well-known and clinically trusted wound care brands.
Start with Cardinal Health.

Table of contents

Antimicrobial dressings (AMD)

AMD

Curity™ AMD Antimicrobial Packing Strips _____	2
Curity™ AMD Antimicrobial Woven Sponges _____	2
Excilon™ AMD Antimicrobial Non-Woven, IV and Drain Sponges _____	4
Kendall™ AMD Antimicrobial Foam Dressings _____	4-6

Advanced wound care

Alginates

Kendall™ Calcium Alginate Dressings _____	8
Cardinal Health™ Reinforced Gelling Fibers (Alginate) _____	8

Composite Dressings

Cardinal Health™ Composite Dressings _____	9
--	---

Contact Layers

Cardinal Health™ Silicone Contact Layers _____	9
Cardinal Health™ Petrolatum Emulsion Contact Layers _____	9

Foam Dressings

Kendall™ AMD Antimicrobial Foam Dressings _____	10
Kendall™ Foam Dressings _____	11
Kendall™ Foam Border Gentle Silicone Adhesion Dressings _____	11
Kendall™ Silicone Foams _____	12
Kendall™ Silicone Bordered and Non-Bordered Foams Lite _____	12
Kendall™ Silicone Bordered Foams _____	12

Hydrocolloids

Kendall™ Alginate Hydrocolloid Dressings _____	13
Cardinal Health™ Hydrocolloids _____	13
Cardinal Health™ Thin Hydrocolloids _____	13

Advanced wound care *(continued)*

Hydrogel Dressings

Kendall™ Amorphous Gel _____	14
Cardinal Health™ Hydrogel _____	14
Kendall™ Hydrogel Impregnated Gauze _____	14
Kendall™ Hydrogel Wound Dressings _____	14

Transparent Dressings

Cardinal Health™ Transparent Film Dressings _____	15
Kendall™ Transparent Film Dressings _____	16

Wound Solution

Kendall™ Sterile Saline Wound Solution _____	16
Cardinal Health™ Wound Cleanser _____	16

Reimbursement information _____ 17-19

Index _____ 20-21



Antimicrobial Dressings with PHMB

Consider Cardinal Health™ AMD Antimicrobial Dressings — with the power of proven, affordable PHMB (polyhexamethylene biguanide) — when your patients need it most. It's protection that works as hard as you do to prevent infection. From Cardinal Health, the pioneer in skin protection.

Hard on bacteria. Easy on skin.

The Cardinal Health line of antimicrobial dressings (AMD) is impregnated with 0.2% PHMB (polyhexamethylene biguanide), a bacteria-killing polymer. So you get powerful antimicrobial action that stays gentle on skin.

- PHMB kills a broad-spectrum of bacteria
- Provides protection against 15 gram-negative, gram-positive and fungi/yeast microorganisms
- Resists bacterial colonization and reduces bacterial penetration through the dressing
- Proven to perform without negatively impacting epithelization
- Low to minimal cytotoxicity
- No known resistances



PHMB
Power to Protect



Curity™ AMD Antimicrobial Packing Strips

Code	Dimensions	Sterility	Pkg. string	Case qty.	Pkg.
7831AMD	¼ in. x 1 yd	Sterile	1 ea/po, 10 po/bx, 5 bx/cs	50	Soft Pouch
7832AMD	½ in. x 1 yd	Sterile	1 ea/po, 10 po/bx, 5 bx/cs	50	Soft Pouch
7833AMD	1 in. x 1 yd	Sterile	1 ea/po, 10 po/bx, 5 bx/cs	50	Soft Pouch

Curity™ AMD Antimicrobial Woven Sponges

Code	Dimensions	Ply	Sterility	Pkg. string	Case qty.	Pkg.
2506-	2 in. x 2 in.	8 ply	Sterile	2 ea/po, 25 po/bx, 30 bx/cs	3,000	Peel Back Package
2533	4 in. x 4 in.	12 ply	Sterile	2 ea/po, 25 po/bx, 24 bx/cs	1,200	Peel Back Package
2539-	4 in. x 4 in.	12 ply	Sterile	10 ea/tr, 128 tr/cs	1,280	Plastic Tray

More power to protect. More power to prevent. More power to heal.

Learn how PHMB works.

PHMB attacks bacteria on and within the dressing to help you keep bacteria out of the wound, limit cross-contamination and promote a healing environment.

Effective against a broad spectrum of bacteria:

Methicillin-resistant Staphylococcus aureus (MRSA)

Vancomycin-resistant Enterococcus (VRE)

Acinetobacter baumannii

Staphylococcus epidermidis

Pseudomonas aeruginosa

Escherichia coli

Candida albicans

Coagulase-negative Staphylococcus

Proteus mirabilis

Serratia marcescens

Enterobacter cloacae

Klebsiella pneumoniae

Enterococcus faecalis

Propionibacterium acnes

Give more patients the prevention they need:

- diabetic patients
- surgical patients
- elderly patients
- CABG and orthopedic surgeries
- caesarean sections
- heel or foot ulcers



PHMB

Power to Protect

Cardinal Health has the brands you trust for protection:

Kerlix™ AMD Rolls, Gauze and Super Sponges

Telfa™ AMD Non-adherent Dressings

Excilon™ AMD Drain and IV Sponges

Kendall™ AMD Foams and Foam Discs

Curity™ AMD Sponges and Packing Strips



Look for this icon to find Cardinal Health AMD products with PHMB throughout this catalog.



Excilon™ AMD Antimicrobial Non-Woven, IV and Drain Sponges

Code	Style	Dimensions	Ply	Sterility	Pkg. string	Case qty.	Pkg.
7089-	IV Sponge	2 in. x 2 in.	6 ply	Sterile	2 ea/po, 35 po/bx, 20 bx/cs	1,400	Peel Back Package
7088	Drain Sponge	4 in. x 4 in.	6 ply	Sterile	2 ea/po, 25 po/bx, 12 bx/cs	600	Peel Back Package



Kendall™ AMD Antimicrobial Foam Dressings

Foam Dressings

Code	Dimensions	Attributes	Sterility	Pkg. string	Case qty.
55522AMD	2 in. x 2 in.		Sterile	1 ea/po, 25 po/bx, 4 bx/cs	100
55535AMD	3 ½ in. x 3 in.	Fenestrated	Sterile	1 ea/po, 10 po/bx, 5 bx/cs	50
55535PAMD	3 ½ in. x 3 in.	Fenestrated with Topsheet	Sterile	1 ea/po, 10 po/bx, 5 bx/cs	50
55544PAMD	4 in. x 4 in.	with Topsheet	Sterile	1 ea/po, 10 po/bx, 5 bx/cs	50
55548AMD	4 in. x 8 in.		Sterile	1 ea/po, 10 po/bx, 5 bx/cs	50
55566AMD	6 in. x 6 in.		Sterile	1 ea/po, 10 po/bx, 5 bx/cs	50
55588AMD	8 in. x 8 in.		Sterile	1 ea/po, 10 po/bx, 5 bx/cs	50
55544AMD	4 in. x 4 in.		Sterile	1 ea/po, 10 po/bx, 5 bx/cs	50

Foam Border Dressings

Code	Dimensions	Attributes	Sterility	Pkg. string	Case qty.
55588BAMD	7 ½ in. x 7 ½ in.	Acrylic Border	Sterile	1 ea/po, 10 po/bx, 5 bx/cs	50
55546BAMD	3 ½ in. x 5 ½ in.	Acrylic Border	Sterile	1 ea/po, 10 po/bx, 5 bx/cs	50
55523BAMD	1 ¾ in. x 3 ¼ in.	Acrylic Border	Sterile	1 ea/po, 10 po/bx, 5 bx/cs	50
55566BAMD	5 ½ in. x 5 ½ in.	Acrylic Border	Sterile	1 ea/po, 10 po/bx, 5 bx/cs	50
55544BAMD	3 ½ in. x 3 ½ in.	Acrylic Border	Sterile	1 ea/po, 10 po/bx, 5 bx/cs	50



Kendall™ AMD Antimicrobial Foam Dressings (continued)

Foam Border Gentle Adhesion Dressings

Code	Dimensions	Attributes	Size	Sterility	Pkg. string	Case qty.
55523BAMDG	2 in. x 3 in.	Silicone Border		Sterile	1 ea/po, 10 po/bx, 5 bx/cs	50
55544BAMDG	4 in. x 4 in.	Silicone Border		Sterile	1 ea/po, 10 po/bx, 5 bx/cs	50
55546BAMDG	4 in. x 6 in.	Silicone Border		Sterile	1 ea/po, 10 po/bx, 5 bx/cs	50
55566BAMDG	6 in. x 6 in.	Silicone Border		Sterile	1 ea/po, 10 po/bx, 5 bx/cs	50
55577BAMDG	7 ½ in. x 7 ½ in.	Silicone Border	Sacrum Small	Sterile	1 ea/po, 10 po/bx, 5 bx/cs	50
55588BAMDG	8 in. x 8 in.	Silicone Border		Sterile	1 ea/po, 10 po/bx, 5 bx/cs	50
55599BAMDG	9 ¼ in. x 9 ¼ in.	Silicone Border	Sacrum Large	Sterile	1 ea/po, 10 po/bx, 5 bx/cs	50

Foam Discs

Code	Dimensions	Sterility	Pkg. string	Case qty.
55511AMD	1 in. Diameter, 4mm Hole	Sterile	1 ea/po, 10 po/bx, 4 bx/cs	40
55512AMD	1 in. Diameter, 7mm Hole	Sterile	1 ea/po, 10 po/bx, 4 bx/cs	40



Kerlix™ AMD Gauze Bandage Rolls

Code	Dimensions	Sterility	Pkg. string	Case qty.	Pkg.
3332	4 ½ in. x 4.1 yds	Sterile	1 rl/po, 100 po/cs	100	Soft Pouch
3331-	4 ½ in. x 4.1 yds	Sterile	1 rl/po, 60 po/cs	60	Rigid Plastic Tray

Kerlix™ AMD Super Sponges

Code	Dimensions	Ply	Sterility	Pkg. string	Case qty.	Pkg.
6662	Medium	6 in. x 6 ¾ in.	Sterile	2 ea/po, 20 po/bx, 12 bx/cs	480	Soft Pouch
6665	Medium	6 in. x 6 ¾ in.	Sterile	5 ea/po, 10 po/bx, 12 bx/cs	600	Soft Pouch
6660	Medium	6 in. x 6 ¾ in.	Sterile	10 ea/tr, 48 tr/cs	480	Plastic Tray

Telfa™ AMD Non-Adherent Dry Dressings

Code	Dimensions	Style	Sterility	Pkg. string	Case qty.
7662-	3 in. x 4 in.	Telfa AMD Pad	Sterile	1 ea/po, 50 po/bx, 18 bx/cs	900
7663	3 in. x 8 in.	Telfa AMD Pad	Sterile	1 ea/po, 50 po/bx, 12 bx/cs	600
7665	4 in. x 5 in.	Telfa AMD Island Dressing	Sterile	1 ea/po, 25 po/bx, 8 bx/cs	200
7666	4 in. x 8 in.	Telfa AMD Island Dressing	Sterile	1 ea/po, 25 po/bx, 4 bx/cs	100
7667	4 in. x 10 in.	Telfa AMD Island Dressing	Sterile	1 ea/po, 25 po/bx, 4 bx/cs	100
7668	4 in. x 14 in.	Telfa AMD Island Dressing	Sterile	1 ea/po, 25 po/bx, 2 bx/cs	50



Advanced wound care



Look for this icon to quickly identify products containing PHMB — your shield against infection

Alginates



Kendall™ Calcium Alginate Dressings

- Dressing absorbs exudate, forms a gel that keeps the wound moist and protects from contamination
- Alginate dressing absorbs wound exudate
- Absorbs approximately 12 times the dressing's weight (in exudate)
- Calcium Alginate composition and weaving pattern aids in maintaining dressing integrity during removal
- Not made with natural rubber latex

Calcium Alginate Dressings

Code	Dimensions	Attributes	Sterility	Pkg. string	Case qty.
9232	2 in. x 2 in.		Sterile	1 ea/po, 10 po/bx, 10 bx/cs	100
9233-	4 in. x 4 in.		Sterile	1 ea/po, 10 po/bx, 5 bx/cs	50
9235	4¼ x 4¼ in.		Sterile	1 ea/po, 100 po/cs	100
9236	4 in. x 4 in.	Plus	Sterile	1 ea/po, 10 po/bx, 5 bx/cs	50
9240	4 in. x 5½ in.		Sterile	1 ea/po, 10 po/bx, 5 bx/cs	50
9238	4 in. x 8 in.		Sterile	1 ea/po, 5 po/bx, 10 bx/cs	50
9239	6 in. x 10 in.		Sterile	1 ea/po, 10 po/bx, 5 bx/cs	50
9242	12 in. x 24 in.		Sterile	1 ea/po, 5 po/cs	5
9231-	12 in.	Rope	Sterile	1 ea/po, 5 po/bx, 4 bx/cs	20
9243	24 in.	Rope	Sterile	1 ea/po, 5 po/bx, 4 bx/cs	20
9244	36 in.	Rope	Sterile	1 ea/po, 5 po/bx, 4 bx/cs	20

Zinc Calcium Alginate Dressings

Code	Dimensions	Attributes	Sterility	Pkg. string	Case qty.
9354-	2 in. x 2 in.	with Zinc	Sterile	1 ea/po, 10 po/bx, 10 bx/cs	100
9355-	4 in. x 4 in.	with Zinc	Sterile	1 ea/po, 10 po/bx, 5 bx/cs	50
9356-	4 in. x 8 in.	with Zinc	Sterile	1 ea/po, 5 po/bx, 10 bx/cs	50
9351-	12 in.	Rope with Zinc	Sterile	1 ea/po, 5 po/bx, 4 bx/cs	20



Cardinal Health™ Reinforced Gelling Fibers (Alginate)

- Non-woven calcium/sodium alginate consisting of a high M alginate in a cohesive fiber construction
- Our reinforced gelling fibers assist in maintaining a moist environment for optimal wound healing and the formation of granulation tissues
- Highly absorbent gelling fibers wick and retain exudate away from the wound bed
- Cohesive fiber construction facilitates easy intact/one piece atraumatic removal
- The silver ions in our reinforced gelling Fibers+Ag protect the dressing from a broad spectrum of microorganisms over a period of up to 14 days
- Not made with natural rubber latex

Reinforced Gelling Fibers

Code	Dimensions	Attributes	Sterility	Pkg. string	Case qty.
RGFB22	2 in. x 2 in.		Sterile	1 ea/po, 10 po/bx, 10 bx/cs	100
RGFB45	4 in. x 4 ¾ in.		Sterile	1 ea/po, 10 po/bx, 10 bx/cs	100
RGFB66	6 in. x 6 in.		Sterile	1 ea/po, 5 po/bx, 10 bx/cs	50
RGFBROPE	18 in.	Rope	Sterile	1 ea/po, 5 po/bx, 10 bx/cs	50

Reinforced Gelling Fibers +AG

Code	Dimensions	Attributes	Sterility	Pkg. string	Case qty.
RGFB22AG	2 in. x 2 in.	+AG	Sterile	1 ea/po, 10 po/bx, 10 bx/cs	100
RGFB45AG	4 in. x 4 ¾ in.	+AG	Sterile	1 ea/po, 10 po/bx, 10 bx/cs	100
RGFB66AG	6 in. x 6 in.	+AG	Sterile	1 ea/po, 5 po/bx, 10 bx/cs	50
RGFB812AG	8 in. x 12 in.	+AG	Sterile	1 ea/po, 5 po/bx, 10 bx/cs	50
RGFBROPEAG	18 in.	Rope +AG	Sterile	1 ea/po, 5 po/bx, 10 bx/cs	50

Composite Dressings



Cardinal Health™ Composite Dressings

- Made with a soft non-woven fabric coated with hypoallergenic acrylic adhesive and an absorbent pad with a low adherent wound contact layer
- Not made with natural rubber latex

Code	Dimensions	Sterility	Pkg. string	Case qty.
COMP23	2 in. x 3 in.	Sterile	1 ea/po, 10 po/bx, 20 bx/cs	200
COMP34	3.2 in. x 4 in.	Sterile	1 ea/po, 10 po/bx, 20 bx/cs	200
COMP36	3.2 in. x 6 in.	Sterile	1 ea/po, 10 po/bx, 20 bx/cs	200
COMP44	4 in. x 4 in.	Sterile	1 ea/po, 10 po/bx, 20 bx/cs	200
COMP48	4 in. x 8 in.	Sterile	1 ea/po, 10 po/bx, 20 bx/cs	200
COMP410	4 in. x 10 in.	Sterile	1 ea/po, 10 po/bx, 20 bx/cs	200
COMP412	4 in. x 12 in.	Sterile	1 ea/po, 10 po/bx, 20 bx/cs	200
COMP414	4 in. x 14 in.	Sterile	1 ea/po, 10 po/bx, 20 bx/cs	200
COMP66	6 in. x 6 in.	Sterile	1 ea/po, 10 po/bx, 20 bx/cs	200

Contact Layers



Cardinal Health™ Silicone Contact Layers

- Wound dressings made with a conformable, open mesh structure and gentle silicone adhesive
- Used to facilitate fluid transfer and protection to the wound bed during dressing changes
- Not made with natural rubber latex

Code	Dimensions	Sterility	Pkg. string	Case qty.
SCL23	2 in. x 3 in.	Sterile	1 ea/po, 10 po/bx, 5 bx/cs	50
SCL34	3 in. x 4 in.	Sterile	1 ea/po, 10 po/bx, 5 bx/cs	50
SCL47	4 in. x 7 in.	Sterile	1 ea/po, 10 po/bx, 5 bx/cs	50
SCL812	8 in. x 12 in.	Sterile	1 ea/po, 5 po/bx, 5 bx/cs	25



Cardinal Health™ Petrolatum Emulsion Contact Layers

- Wound dressings made with a flexible, non-adherent knitted cellulose acetate fabric coated with a formulated petrolatum emulsion
- Used to facilitate fluid transfer and protection to the wound bed during dressing changes
- Not made with natural rubber latex

Code	Dimensions	Sterility	Pkg. string	Case qty.
PETROCL33	3 in. x 3 in.	Sterile	1 ea/po, 50 po/bx, 12 bx/cs	600
PETROCL38	3 in. x 8 in.	Sterile	1 ea/po, 24 po/bx, 6 bx/cs	144

Foam Dressings



Kendall™ AMD Antimicrobial Foam Dressings

- Impregnated with PHMB (polyhexamethylene biguanide)
- Kills a broad-spectrum of bacteria
- Provides protection against 15 gram-negative, gram-positive and fungi/yeast microorganisms
- Resists bacterial colonization and reduces bacterial penetration through the dressing
- Proven to perform without negatively impacting epithelization
- Low to minimal cytotoxicity
- No known resistances
- Hypoallergenic
- Not made with natural rubber latex



Antimicrobial Foam Dressings

Code	Dimensions	Attributes	Sterility	Pkg. string	Case qty.
55522AMD	2 in. x 2 in.		Sterile	1 ea/po, 25 po/bx, 4 bx/cs	100
55535AMD	3 ½ in. x 3 in.	Fenestrated	Sterile	1 ea/po, 10 po/bx, 5 bx/cs	50
55535PAMD	3 ½ in. x 3 in.	Fenestrated with Topsheet	Sterile	1 ea/po, 10 po/bx, 5 bx/cs	50
55544PAMD	4 in. x 4 in.	with Topsheet	Sterile	1 ea/po, 10 po/bx, 5 bx/cs	50
55544AMD	4 in. x 4 in.		Sterile	1 ea/po, 10 po/bx, 5 bx/cs	50
55548AMD	4 in. x 8 in.		Sterile	1 ea/po, 10 po/bx, 5 bx/cs	50
55566AMD	6 in. x 6 in.		Sterile	1 ea/po, 10 po/bx, 5 bx/cs	50
55588AMD	8 in. x 8 in.		Sterile	1 ea/po, 10 po/bx, 5 bx/cs	50

Antimicrobial Foam Border Dressings

Code	Dimensions	Attributes	Sterility	Pkg. string	Case qty.
55523BAMD	1 ¾ in. x 3 ¼ in.	Acrylic Adhesive	Sterile	1 ea/po, 10 po/bx, 5 bx/cs	50
55544BAMD	3 ½ in. x 3 ½ in.	Acrylic Adhesive	Sterile	1 ea/po, 10 po/bx, 5 bx/cs	50
55546BAMD	3 ½ in. x 5 ½ in.	Acrylic Adhesive	Sterile	1 ea/po, 10 po/bx, 5 bx/cs	50
55566BAMD	5 ½ in. x 5 ½ in.	Acrylic Adhesive	Sterile	1 ea/po, 10 po/bx, 5 bx/cs	50
55588BAMD	7 ½ in. x 7 ½ in.	Acrylic Adhesive	Sterile	1 ea/po, 10 po/bx, 5 bx/cs	50

Antimicrobial Foam Border Gentle Adhesion Dressings

Code	Dimensions	Attributes	Size	Sterility	Pkg. string	Case qty.
55523BAMDG	2 in. x 3 in.	Silicone Adhesive		Sterile	1 ea/po, 10 po/bx, 5 bx/cs	50
55544BAMDG	4 in. x 4 in.	Silicone Adhesive		Sterile	1 ea/po, 10 po/bx, 5 bx/cs	50
55546BAMDG	4 in. x 6 in.	Silicone Adhesive		Sterile	1 ea/po, 10 po/bx, 5 bx/cs	50
55566BAMDG	6 in. x 6 in.	Silicone Adhesive		Sterile	1 ea/po, 10 po/bx, 5 bx/cs	50
55577BAMDG	7 ½ in. x 7 ½ in.	Silicone Adhesive	Sacrum Small	Sterile	1 ea/po, 10 po/bx, 5 bx/cs	50
55588BAMDG	8 in. x 8 in.	Silicone Adhesive		Sterile	1 ea/po, 10 po/bx, 5 bx/cs	50
55599BAMDG	9 ¼ in. x 9 ¼ in.	Silicone Adhesive	Sacrum Large	Sterile	1 ea/po, 10 po/bx, 5 bx/cs	50

Antimicrobial Foam Discs

Code	Dimensions	Sterility	Pkg. string	Case qty.
55511AMD	1 in. Diameter, 4mm Hole	Sterile	1 ea/po, 10 po/bx, 4 bx/cs	40
55512AMD	1 in. Diameter, 7mm Hole	Sterile	1 ea/po, 10 po/bx, 4 bx/cs	40

Foam Dressings



Kendall™ Foam Dressings

- Soft hydrophilic foam protects and cushions wounds
- Top-sheet material found on Kendall™ foam dressing with topsheet and island dressing allows moisture vapor transmission while preventing strike through
- Not made with natural rubber latex

Foam Dressings

Code	Dimensions	Attributes	Sterility	Pkg. string	Case qty.
55522	2 in. x 2 in.		Sterile	1 ea/po, 25 po/bx, 4 bx/cs	100
55533	3 in. x 3 in.		Sterile	1 ea/po, 10 po/bx, 5 bx/cs	50
55535	3 ½ in. x 3 in.	Fenestrated	Sterile	1 ea/po, 10 po/bx, 5 bx/cs	50
55544	4 in. x 4 in.		Sterile	1 ea/po, 10 po/bx, 5 bx/cs	50
55548	4 in. x 8 in.		Sterile	1 ea/po, 10 po/bx, 5 bx/cs	50
55555	5 in. x 5 in.		Sterile	1 ea/po, 10 po/bx, 5 bx/cs	50
55566	6 in. x 6 in.		Sterile	1 ea/po, 10 po/bx, 5 bx/cs	50
55588	8 in. x 8 in.		Sterile	1 ea/po, 10 po/bx, 5 bx/cs	50

Foam Dressings with Topsheet

Code	Dimensions	Attributes	Sterility	Pkg. string	Case qty.
55522P	2 in. x 2 in.	with Topsheet	Sterile	1 ea/po, 25 po/bx, 4 bx/cs	100
55533P	3 in. x 3 in.	with Topsheet	Sterile	1 ea/po, 10 po/bx, 5 bx/cs	50
55535P	3 ½ in. x 3 in.	Fenestrated with Topsheet	Sterile	1 ea/po, 10 po/bx, 5 bx/cs	50
55544P	4 in. x 4 in.	with Topsheet	Sterile	1 ea/po, 10 po/bx, 5 bx/cs	50
55548P	4 in. x 8 in.	with Topsheet	Sterile	1 ea/po, 10 po/bx, 5 bx/cs	50
55566P	6 in. x 6 in.	with Topsheet	Sterile	1 ea/po, 10 po/bx, 5 bx/cs	50
55588P	8 in. x 8 in.	with Topsheet	Sterile	1 ea/po, 10 po/bx, 5 bx/cs	50

Foam Island Dressings

Code	Dimensions	Attributes	Sterility	Pkg. string	Case qty.
55544B	4 in. x 4 in. Dressing, 2 in. x 2 in. Pad	Acrylic Adhesive Border	Sterile	1 ea/po, 10 po/bx, 5 bx/cs	50
55566B	6 in. x 6 in. Dressing, 4 in. x 4 in. Pad	Acrylic Adhesive Border	Sterile	1 ea/po, 10 po/bx, 5 bx/cs	50
55588B	8 in. x 8 in. Dressing, 6 in. x 6 in. Pad	Acrylic Adhesive Border	Sterile	1 ea/po, 10 po/bx, 5 bx/cs	50



Kendall™ Foam Border Gentle Silicone Adhesion Dressings

- Gentle silicone adhesive minimizes skin trauma and remains adherent after adjusting the dressing
- High MVTR (Moisture Vapor Transmission Rate) backing reduces friction
- Small pore size helps to prevent microbe and cell in-growth, while also providing a soft feel and to help relieve pressure (*Data on file)
- Not made with natural rubber latex

Code	Dimensions	Attributes	Size	Sterility	Pkg. string	Case qty.
55523BG	1 ¾ in. x 3 ¼ in.	Silicon Border		Sterile	1 ea/po, 10 po/bx, 5 bx/cs	50
55544BG	3 ½ in. x 3 ½ in.	Silicon Border		Sterile	1 ea/po, 10 po/bx, 5 bx/cs	50
55546BG	3 ½ in. x 5 ½ in.	Silicon Border		Sterile	1 ea/po, 10 po/bx, 5 bx/cs	50
55566BG	5 ½ in. x 5 ½ in.	Silicon Border		Sterile	1 ea/po, 10 po/bx, 5 bx/cs	50
55588BG	7 ½ in. x 7 ½ in.	Silicon Border		Sterile	1 ea/po, 10 po/bx, 5 bx/cs	50
55577BG	7 ½ in. x 7 ½ in.	Silicon Border	Sacrum Small	Sterile	1 ea/po, 10 po/bx, 5 bx/cs	50
55599BG	9 ¼ in. x 9 ¼ in.	Silicon Border	Sacrum Large	Sterile	1 ea/po, 10 po/bx, 5 bx/cs	50

Foam Dressings



Kendall™ Silicone Foams

- Highly absorbent foam pad, protective film backing and gentle silicone adhesive across the face of the dressing
- Protects and absorbs exudate from the wound
- Helps facilitate a moist environment to promote healing
- Not made with natural rubber latex

Code	Dimensions	Attributes	Sterility	Pkg. string	Case qty.
FMTRACH	3 ½ in. x 3 ½ in.	Non-Adhesive Foam Tracheostomy	Sterile	1 ea/po, 5 po/bx, 5 bx/cs	25
FM44	4 in. x 4 in.		Sterile	1 ea/po, 10 po/bx, 5 bx/cs	50
FM48	4 in. x 8 in.		Sterile	1 ea/po, 10 po/bx, 5 bx/cs	50
FMHEEL	4.7 in. x 7 ½ in.	Heel	Sterile	1 ea/po, 5 po/bx, 5 bx/cs	25
FM66	6 in. x 6 in.		Sterile	1 ea/po, 10 po/bx, 5 bx/cs	50
FM88	8 in. x 8 in.		Sterile	1 ea/po, 5 po/bx, 5 bx/cs	25

Kendall™ Silicone Bordered and Non-Bordered Foams Lite

- Foam dressings made within absorbent thin foam pad, protective film backing and gentle silicone adhesive
- Used to absorb and retain exudate away from the wound and facilitate a moist environment
- Gentle silicone adhesive helps allow for atraumatic removal
- Not made with natural rubber latex

Code	Dimensions	Attributes	Sterility	Pkg. string	Case qty.
BFM22LTE	1.6 in. x 2 in.	Bordered	Sterile	1 ea/po, 10 po/bx, 5 bx/cs	50
BFM25LTE	2 in. x 5 in.	Bordered	Sterile	1 ea/po, 10 po/bx, 5 bx/cs	50
BFM33LTE	3 in. x 3 in.	Bordered	Sterile	1 ea/po, 10 po/bx, 5 bx/cs	50
BFM44LTE	4 in. x 4 in.	Bordered	Sterile	1 ea/po, 10 po/bx, 5 bx/cs	50
BFM66LTE	6 in. x 6 in.	Bordered	Sterile	1 ea/po, 10 po/bx, 5 bx/cs	50
FM23LTE	2.4 in. x 3.4 in.	Non-Bordered	Sterile	1 ea/po, 10 po/bx, 5 bx/cs	50
FM44LTE	4 in. x 4 in.	Non-Bordered	Sterile	1 ea/po, 10 po/bx, 5 bx/cs	50

Kendall™ Silicone Bordered Foams

- Multi-layered bordered foam dressings consist of a highly absorbent foam pad, additional layer containing super absorbent polyacrylate fibers, protective film backing, gentle silicone adhesive and a secure border
- Used to absorb and retain exudate away from the wound and facilitate a moist environment
- Gentle silicone adhesive helps allow for atraumatic removal
- Not made with natural rubber latex

Code	Dimensions	Size	Sterility	Pkg. string	Case qty.
BFM33	3 in. x 3 in.		Sterile	1 ea/po, 10 po/bx, 5 bx/cs	50
BFM44	4 in. x 4 in.		Sterile	1 ea/po, 10 po/bx, 5 bx/cs	50
BFM48	4 in. x 8 in.		Sterile	1 ea/po, 5 po/bx, 5 bx/cs	25
BFM412	4 in. x 12 in.		Sterile	1 ea/po, 5 po/bx, 5 bx/cs	25
BFM66	6 in. x 6 in.		Sterile	1 ea/po, 10 po/bx, 5 bx/cs	50
BFM68	6 in. x 8 in.		Sterile	1 ea/po, 5 po/bx, 5 bx/cs	25
BFMSMSCRL	7.2 in. x 7.2 in.	Sacrum Small	Sterile	1 ea/po, 5 po/bx, 5 bx/cs	25
BFMLGSCRL	9.2 in. x 9.2 in.	Sacrum Large	Sterile	1 ea/po, 5 po/bx, 5 bx/cs	25

Hydrocolloids



Kendall™ Alginate Hydrocolloid Dressings

- Wound dressing designed for use on light to moderately exuding wounds
- Provides an optimal environment for moist wound healing
- Alginate formulation provides greater absorbency resulting in fewer dressing changes
- Thin profile aids in product conformability
- Not made with natural rubber latex

Code	Dimensions	Attributes	Sterility	Pkg. string	Case qty.
9807	2 ½ in. x 2 ½ in.	Border	Sterile	1 ea/po, 5 po/bx, 10 bx/cs	50
9808	4 in. x 4 in.	Border	Sterile	1 ea/po, 5 po/bx, 10 bx/cs	50
9801-	4 in. x 4 in.		Sterile	1 ea/po, 5 po/bx, 10 bx/cs	50
9805	4 in. x 5 in.	Sacral	Sterile	1 ea/po, 5 po/bx, 10 bx/cs	50
9802	6 in. x 6 in.		Sterile	1 ea/po, 5 po/bx, 10 bx/cs	50
9809-	6 in. x 6 in.	Border	Sterile	1 ea/po, 5 po/bx, 10 bx/cs	50
9806	6 in. x 7 in.	Sacral	Sterile	1 ea/po, 5 po/bx, 10 bx/cs	50
9804	8 in. x 8 in.		Sterile	1 ea/po, 5 po/bx, 10 bx/cs	50



Cardinal Health™ Hydrocolloids

- Helps facilitate local wound healing by absorbing wound exudate and creating a moist environment
- Proprietary gelling formulation allows for most of the gel to be removed together with the dressing, resulting in little or no damage to the newly formed tissue
- Flexible, low-friction film backing reduces skin friction with clothing and patient bedding
- Moisture-resistant and will remain in place while bathing
- Not made with natural rubber latex

Code	Dimensions	Attributes	Sterility	Pkg. string	Case qty.
HC44	4 in. x 4 in.		Sterile	1 ea/po, 10 po/bx, 20 bx/cs	200
HCOVAL	4 in. x 4 ¾ in.	Oval	Sterile	1 ea/po, 10 po/bx, 20 bx/cs	200
HC66	6 in. x 6 in.		Sterile	1 ea/po, 5 po/bx, 16 bx/cs	80
HCTRI	6 in. x 7 in.	Triangle	Sterile	1 ea/po, 5 po/bx, 16 bx/cs	80
HC88	8 in. x 8 in.		Sterile	1 ea/po, 5 po/bx, 16 bx/cs	80



Cardinal Health™ Thin Hydrocolloids

- Thin profile design with rounded corners to help reduce dressing roll-up during product use
- Semi-translucent film backing allows for easier monitoring of the healing process without dressing removal
- Not made with natural rubber latex

Code	Dimensions	Attributes	Sterility	Pkg. string	Case qty.
HCTSPOTS	1 ¾ in. x 1 ½ in.	Spots	Sterile	1 ea/po, 10 po/bx, 20 bx/cs	200
HCT33	3 in. x 3 in.		Sterile	1 ea/po, 10 po/bx, 20 bx/cs	200
HCT44	4 in. x 4 in.		Sterile	1 ea/po, 10 po/bx, 20 bx/cs	200
HCTOVAL	4 in. x 4 ¾ in.	Oval	Sterile	1 ea/po, 10 po/bx, 20 bx/cs	200
HCT66	6 in. x 6 in.		Sterile	1 ea/po, 5 po/bx, 16 bx/cs	80

Hydrogel Dressings



Kendall™ Amorphous Gel

- Clear wound gel dressing helps provide the optimal environment for natural wound healing
- Formulated for the management of Stage II-IV pressure ulcers, stasis ulcers, first and second degree burns, cuts, abrasions and minor irritations of the skin
- Glycerin-based hydrogel will not dry out
- Not made with natural rubber latex

Code	Size	Sterility	Pkg. string	Case qty.
9250-	½ oz.	Non-Sterile	1 ea/bx, 24 bx/cs	24
9251-	1 oz.	Non-Sterile	1 ea/bx, 12 bx/cs	12
9252	3 oz.	Non-Sterile	1 ea/bx, 12 bx/cs	12



Cardinal Health™ Hydrogel

- Clean, greaseless hydrogel used for the maintenance of a moist wound environment
- Excellent viscosity helps aid product application
- Helps maintain a moist wound environment
- Hydrogel+Ag contains an antimicrobial silver compound which is an effective barrier to bacterial penetration that inhibits the growth of a broad spectrum of microorganisms that come into contact with the gel
- The hydrogel+Ag formulation is effective in assisting the debridement and desloughing process in dry necrotic wounds
- Not made with natural rubber latex

Hydrogel + AG

Code	Size	Attributes	Sterility	Pkg. string	Case qty.
HG15OZAG	1.5 oz	+AG	Non-Sterile	1 ea/bx, 10 bx/ct, 4 ct/cs	40



Kendall™ Hydrogel Impregnated Gauze

- Clear hydrogel impregnated in a gauze sponge
- Provides optimal healing environment for moist wound healing
- Extends dressing wear times
- Not made with natural rubber latex

Code	Dimensions	Sterility	Pkg. string	Case qty.
9255	2 in. x 2 in.	Sterile	1 ea/po, 25 po/bx, 4 bx/cs	100
9256	4 in. x 4 in.	Sterile	1 ea/po, 25 po/bx, 4 bx/cs	100
9257	8 in. x 4 in.	Sterile	1 ea/po, 25 po/bx, 4 bx/cs	100



Kendall™ Hydrogel Wound Dressings

- Unique disc shape maximizes wound coverage and helps fill shallow cavities
- Protects wounds as it provides a moist environment
- Translucent gel allows wound visualization
- Cooling sensation helps reduce pain upon dressing change
- Non-traumatic removal when changing the dressings
- Not made with natural rubber latex

Code	Dimensions	Sterility	Pkg. string	Case qty.
8884476139	3 in.	Sterile	1 ea/po, 5 po/bx, 10 bx/cs	50
8884476154	4 ¾ in.	Sterile	1 ea/po, 5 po/bx, 10 bx/cs	50

Transparent Dressings



Cardinal Health™ Transparent Film Dressings

- Conform and flex with patient movements
- Moisture vapor permeability allows excess moisture out while maintaining a moist wound healing environment
- Transparency allows observation without removing the dressing
- Acrylic adhesive ensures securement
- Advanced securement with notches, tape and border, reducing tension on stitches, wounds and IV site
- Can be worn up to a week
- Waterproof
- Not made with natural rubber latex

Frame Style Window Out

Code	Dimensions	Sterility	Pkg. string	Case qty.
TD-22C	1.75 in. x 1.75 in.	Sterile	1 ea/po, 100 po/bx, 4 bx/cs	400
TD-24C	2.375 in. x 2.75 in.	Sterile	1 ea/po, 100 po/bx, 4 bx/cs	400
TD-26C	4 in. x 4.75 in.	Sterile	1 ea/po, 50 po/bx, 4 bx/cs	200
TD-30	4 in. x 4.5 in.	Sterile	1 ea/po, 50 po/bx, 4 bx/cs	200

Frame Style Window In

Code	Dimensions	Sterility	Pkg. string	Case qty.
TD-27	4 in. x 10 in.	Sterile	1 ea/po, 20 po/bx, 4 bx/cs	80
TD-28	6 in. x 8 in.	Sterile	1 ea/po, 10 po/bx, 8 bx/cs	80
TD-29	8 in. x 12 in.	Sterile	1 ea/po, 10 po/bx, 8 bx/cs	80

First Aid Style

Code	Dimensions	Sterility	Pkg. string	Case qty.
TD-20	2.375 in. x 2.75 in.	Sterile	1 ea/po, 100 po/bx, 6 bx/cs	600
TD-21	4 in. x 4.75 in.	Sterile	1 ea/po, 50 po/bx, 4 bx/cs	200

IV Style with Border

Code	Dimensions	Sterility	Pkg. string	Case qty.
TD-14	2.375 in. x 2.75 in.	Sterile	1 ea/po, 100 po/bx, 4 bx/cs	400
TD-16	4 in. x 4.75 in.	Sterile	1 ea/po, 50 po/bx, 4 bx/cs	200

IV Style

Code	Dimensions	Sterility	Pkg. string	Case qty.
TD-10	2 in. x 2.25 in.	Sterile	1 ea/po, 100 po/bx, 4 bx/cs	400
TD-33	2.75 in. x 3.25 in.	Sterile	1 ea/po, 100 po/bx, 4 bx/cs	400
TD-35	3.25 in. x 4.25 in.	Sterile	1 ea/po, 50 po/bx, 4 bx/cs	200
TD-55	3.5 in. x 4.5 in.	Sterile	1 ea/po, 50 po/bx, 4 bx/cs	200
TD-50	4 in. x 6.125 in.	Sterile	1 ea/po, 25 po/bx, 4 bx/cs	100

Transparent Dressings



Kendall™ Transparent Film Dressings

- Dressings used for IV sites, donor sites, superficial skin ulcers and minor burns
- Simple, effective delivery system allows one-handed application
- Permeable to oxygen and moisture vapor
- Barrier to bacteria and fluid contaminants
- Not made with natural rubber latex
- Window-style delivery method

Transparent Film Dressings

Code	Dimensions	Sterility	Pkg. string	Case qty.
6651-	1 ½ in. x 1 ½ in.	Sterile	1 ea/po, 100 po/bx, 4 bx/cs	400
6640	2 in. x 2 ¾ in.	Sterile	1 ea/po, 100 po/bx, 4 bx/cs	400
6641	4 in. x 4 ¾ in.	Sterile	1 ea/po, 50 po/bx, 4 bx/cs	200
6647-	4 in. x 8 in.	Sterile	1 ea/po, 20 po/bx, 6 bx/cs	120
6642	6 in. x 8 in.	Sterile	1 ea/po, 10 po/bx, 4 bx/cs	40
6648	8 in. x 10 in.	Sterile	1 ea/po, 10 po/bx, 8 bx/cs	80

Wound Solution



Kendall™ Sterile Saline Wound Solution

- USP normal saline
- Allows for multiple, single patient uses
- Not made with natural rubber latex

Code	Size	Sterility	Pkg. string	Case qty.
90SAL	3 fl oz. (90mL)	Sterile	12 ea/cs	12
210SAL	7.1 fl oz. (210mL)	Sterile	12 ea/cs	12



Cardinal Health™ Wound Cleanser

- Gentle, non-irritating formulation aids in the removal of foreign materials such as dirt and debris from wounds and skin
- No-rinse formula
- Non-ionic surfactant
- Two option trigger sprayer (spray or stream)
- Delivers a PSI of 9.0 at 3-inch spraying distance (which falls within AHCPH guidelines of 4 to 15 PSI)
- Not made with natural rubber latex

Code	Size	Sterility	Pkg. string	Case qty.
WC8OZ	8 oz.	Sterile	12 ea/cs	12

Reimbursement information

If a product code is not listed, a HCPCS code has not been identified.

Item code	Suggested HCPCS*	Item code	Suggested HCPCS*	Item code	Suggested HCPCS*
1109	A6206	2419T	A4452	4032	A6216
1050-	A6402	2438-	A6216	4095-	A9270
1111	A6207	2502-	A9270	4152	A9270
1113	A6208	2506-	A6222	4221-	A9270
1114-	A6208	2531C	A4450	441000	A6402
1115	A6208	2533	A6222	441001	A6402
1167	A6217	2539-	A6222	441002	A6402
1169	A6403	2554	A9270	441009	A6402
1238-	A6403	2556-	A6216	441100	A6445
12388	A6217	2562	A6203	441101	A6446
1272	A6218	2563	A6204	441102	A6443
1418	A9270	2564	A6205	441103	A6446
1596C	A4452	2565	A6205	441105	A6446
1596T	A4452	2585	A6403	441106	A6446
1700	A6216	2634-	A6216	441107	A6446
1713-	A6216	2666	A9270	441108	A6445
1792	A6216	2733-	A6216	441121	A6444
1801-	A6442	2754-	A9270	441122	A6443
1806	A6402	2835-	A6216	441123	A6443
1892	A6443	2847	A9270	441124	A6442
1914C	A4452	2891	A6217	441200	A6403
1914S1	A4452	2913	A6402	441201	A6403
1914T	A4452	2944-	A9270	441202	A6403
1961	A6402	30066	A6218	441203	A6216
2059-	A9270	3027C	A4450	441204	A6402
210SAL	A6260	3033	A6402	441205	A6216
2132-	A6402	3041	A6216	441206	A6216
2146-	A6216	3063C	A4452	441208	A6402
2187-	A6402	3142C	A4452	441209	A6216
2230-	A6448	3157	A6402	441210	A6216
2231-	A6448	3175-	A9270	441211	A6402
2232-	A6449	3208	A6216	441214	A6216
2236	A6449	3267C	A4452	441215	A6216
2238-	A6450	3279	A6218	441216	A6216
2239	A6448	3324	A6443	441217	A6216
2242-	A6448	3331-	A6266	441218	A6216
2244-	A6449	3332	A6266	441219	A6216
2247	A6449	3335-	A6407	441250	A6442
2249	A6450	3354C	A4452	441251	A6443
2252	A6216	3381	A6216	441400	A6216
2259-	A6403	3394C	A4452	441401	A6216
2283-	A9270	3489	A9270	441402	A6216
2304C	A4450	3615C	A4450	441403	A6402
2346-	A6216	3967	A6402	441405	A6402
2394-	A9270	3968--	A6402	441406	A6402
2419C	A4452	397110	A6402	441407	A6402
2419S2	A4452	397310	A6402	441408	A6402

Item code	Suggested HCPCS*
441412	A6402
441500	A6448
441501	A6449
441502	A6449
441503	A6450
441504	A6448
441505	A6449
441506	A6449
441507	A6450
441601	A6216
441602	A6403
441603	A6217
441604	A6403
441605	A6217
441606	A6403
441607	A6217
442211	A6216
442212	A6402
442213	A6402
442214	A6402
442215	A6402
442308	A6216
4519-	A9270
4588	A6402
5042-	A6216
5072	A6216
55511AMD	A4221 OR A4649
55512AMD	A4221 OR A4649
55522	A6209
55522AMD	A6209
55522P	A6209
55523BAMD	A6212
55523BG	A6209
55533	A6209
55533P	A6209
55535	A6209
55535AMD	A6209
55535P	A6209
55535PAMD	A6209
55544	A6209
55544AMD	A6209
55544B	A6212
55544BAMD	A6212
55544BG	A6212
55544P	A6209
55544PAMD	A6209
55546BAMD	A6212

Item code	Suggested HCPCS*
55546BG	A6212
55548	A6210
55548AMD	A6210
55548P	A6210
55555	A6210
55566	A6210
55566AMD	A6210
55566B	A6212
55566BAMD	A6212
55566BG	A6212
55566P	A6210
55577BAMDG	A4649
55577BG	A4649
55588	A6211
55588AMD	A6211
55588B	A6213
55588BAMD	A6213
55588BG	A6213
55588P	A6211
55599BAMDG	A4649
55599BG	A4649
5806C	A4450
5807C-	A4450
5808C	A4450
6017	A6402
6021	A6442
6035-	A6218
6112	A6222
6113-	A6223
6114-	A6223
6115-	A6222
6116	A6223
6120-	A6402
6121-	A6224
6132	A6216
6196D	A6217
6197D	A6218
6198D	A6218
6241	A6402
6242P	A6402
6244-	A9270
6246-	A9270
6309	A6216
6318-	A6402
6613C	A4450
6640	A6257
6640WF	A6257

Item code	Suggested HCPCS*
6641	A6258
6641WF	A6258
6642	A6258
6642WF	A6258
6647-	A6258
6647WF	A6258
6648	A6259
6648WF	A6259
6651-	A6257
6651WF	A6257
6660	A6223
6662	A6223
6665	A6223
6715-	A6446
6716	A6446
6720-	A6445
6725	A6446
6730	A6446
6735	A6443
68659	A6218
69043	A6218
6939	A6402
7046C	A4450
7083-	A6402
7084-	A6402
7084AMD	A6222
7085	A6402
7086-	A6402
7087-	A6402
7088	A6222
7089-	A6222
7137C	A4452
7138C	A4452
7139C	A4452
7140C	A4452
7196D	A6403
7197D	A6403
7198D	A6404
7267	A6216
7310-	A6403
7540-	A6220
7541	A6220
7542	A6220
7544	A6220
7550-	A6219
7551	A6219
7605-	A6402

Reimbursement information

If a product code is not listed, a HCPCS code has not been identified.

Item code	Suggested HCPCS*
7631	A6407
7632	A6407
7633	A6407
7634	A6407
7643	A6402
7662-	A6222
7663	A6223
7665	A6203
7666	A6203
7667	A6204
7668	A6204
7831	A6266
7831AMD	A6266
7832	A6266
7832AMD	A6266
7833	A6266
7833AMD	A6266
7834	A6266
8042	A6402
8043	A6402
8044--	A6402
8045	A6402
8046-	A6402
8047-	A6402
8190A	A6217
8192A	A6218
8194A	A6218
8533C	A4452
8534C	A4452
8535C	A4452
8536C	A4452
8884431302	A6266
8884431605	A6266
8884432000	A6266
8884433301	A6266
8884433400	A6222
8884433500	A6222
8884433605	A6266
8884436400	A6222
8884437100	A6223
8884437201	A6222
8884437305	A6266
8884476139	A6242
8884476154	A6243
8886834000	A6206

Item code	Suggested HCPCS*
8886834100	A6207
8886834200	A6208
9022	A6216
9023	A6216
9024-	A6216
9026	A6216
905AL	A6260
9132	A6216
9134	A6216
9190A	A6252
9192A	A6253
9194A	A6253
9231-	A6199
9232	A6196
9233-	A6196
9236	A6196
9238	A6197
9239	A6198
9240	A6197
9242	A6198
9243	A6199
9244	A6199
9250-	A6261
9251-	A6261
9252	A6261
9255	A6231
9256	A6231
9257	A6232
9351-	A6199
9354-	A6196
9355-	A6196
9356-	A6197
9411C	A4452
9412C	A4452
9413C	A4452
9414C	A4452
9416C	A4452
9640-	A6402
9642-	A6402
9801-	A6234
9802	A6235
9804	A6236
9805	A6238
9806	A6238
9807	A6237

Item code	Suggested HCPCS*
9808	A6237
9809-	A6238

Index

Code	Page	Code	Page	Code	Page
#					
90SAL	16	9805	13	55599BG	11
210SAL	16	9806	13	8884476139	14
2506-	2	9807	13	8884476154	14
2533	2	9808	13		
2539-	2	9809-	13	B	
3331-	6	55511AMD	6, 10	BFM22LTE	12
3332	6	55512AMD	6, 10	BFM25LTE	12
6640	16	55522	11	BFM33	12
6641	16	55522AMD	4, 10	BFM33LTE	12
6642	16	55522P	11	BFM44	12
6647-	16	55523BAMD	4, 10	BFM44LTE	12
6648	16	55523BAMDG	5, 10	BFM48	12
6651-	16	55523BG	11	BFM66	12
6660	6	55533	11	BFM66LTE	12
6662	6	55533P	11	BFM68	12
6665	6	55535	11	BFM412	12
7088	4	55535AMD	4, 10	BFMLGSCRL	12
7089-	4	55535P	11	BFMSMSCRL	12
7662-	6	55535PAMD	4, 10		
7663	6	55544	11	C	
7665	6	55544AMD	4, 10	COMP23	9
7666	6	55544B	11	COMP34	9
7667	6	55544BAMD	4, 10	COMP36	9
7668	6	55544BAMDG	5, 10	COMP44	9
7831AMD	2	55544BG	11	COMP48	9
7832AMD	2	55544P	11	COMP66	9
7833AMD	2	55544PAMD	4, 10	COMP410	9
9231-	8	55546BAMD	4, 10	COMP412	9
9232	8	55546BAMDG	5, 10	COMP414	9
9233-	8	55546BG	11		
9235	8	55548	11	F	
9236	8	55548AMD	4, 10	FM23LTE	12
9238	8	55548P	11	FM44	12
9239	8	55555	11	FM44LTE	12
9240	8	55566	11	FM48	12
9242	8	55566AMD	4, 10	FM66	12
9243	8	55566B	11	FM88	12
9244	8	55566BAMD	4, 10	FMHEEL	12
9250-	14	55566BAMDG	5, 10	FMTRACH	12
9251-	14	55566BG	11		
9252	14	55566P	11	H	
9255	14	55577BAMDG	5, 10	HC44	13
9256	14	55577BG	11	HC66	13
9257	14	55588	11	HC88	13
9351-	8	55588AMD	4, 10	HCOVAL	13
9354-	8	55588B	11	HCT33	13
9355-	8	55588BAMD	4, 10	HCT44	13
9356-	8	55588BAMDG	5, 10	HCT66	13
9801-	13	55588BG	11	HCTOVAL	13
9802	13	55588P	11		
9804	13	55599BAMDG	5, 10		

Index

Code	Page
HCTRI _____	13
HCTSPOTS _____	13
HG15OZAG _____	14

P

PETROCL33 _____	9
PETROCL38 _____	9

R

RGFB22 _____	8
RGFB22AG _____	8
RGFB45 _____	8
RGFB45AG _____	8
RGFB66 _____	8
RGFB66AG _____	8
RGFB812AG _____	8
RGFBROPE _____	8
RGFBROPEAG _____	8

S

SCL23 _____	9
SCL34 _____	9
SCL47 _____	9
SCL812 _____	9

T

TD-10 _____	15
TD-14 _____	15
TD-16 _____	15
TD-20 _____	15
TD-21 _____	15
TD-22C _____	15
TD-24C _____	15
TD-26C _____	15
TD-27 _____	15
TD-28 _____	15
TD-29 _____	15
TD-30 _____	15
TD-33 _____	15
TD-35 _____	15
TD-50 _____	15
TD-55 _____	15

W

WC8OZ _____	16
-------------	----

