

# Cardinal Health™ Decontamination Glove

## Cardinal Health™ Decontamination

Powder-Free Nitrile Examination Gloves



Nitrile Decontamination exam gloves feature a 16 inch cuff to provide added coverage and protection. Nitrile Decontamination is the thickest glove in our Cardinal Health™ synthetic exam glove portfolio, designed to offer protection for the Central Sterile, Decontamination space where there is potential exposure to sharps and chemicals.

### Features and Benefits:

- Thickest Cardinal Health™ nitrile exam glove in portfolio
- 16 inch extended length
- Full textured fingers
- Tested for use with various lab chemicals<sup>1</sup> and chemotherapy drugs<sup>2</sup>  
*Note: Test results are printed directly on the packaging.*
- Tested for Puncture Resistance<sup>3</sup>



Cat. No.	Size	Pkg.
88NDS	Small	50/bx, 10 bx/cs
88NDM	Medium	50/bx, 10 bx/cs
88NDL	Large	50/bx, 10 bx/cs
88NDXL	X-Large	50/bx, 10 bx/cs

- Powder-free
- Not made with natural rubber latex
- Color: Blue
- Nominal fingertip thickness: 10.9mil

For more information or to order, contact your Cardinal Health sales representative, call 877.CARDINAL or visit [cardinalhealth.com](http://cardinalhealth.com).

[cardinalhealth.com/gloves](http://cardinalhealth.com/gloves)

1. Caution: Gloves used for protection against laboratory chemical exposure should be selected specifically for the type of chemicals being used. Note: Review Safety Data Sheets (SDS) for the chemicals/drugs being used to determine required level of protection.

2. This glove has been tested for resistance to permeation of various chemotherapy drugs per ASTM D6978, "Standard Practice for Assessment of Resistance of Medical Gloves to Permeation by Chemotherapy Drugs."

3. This glove has been tested for puncture resistance per ASTM F1342-91, "Standard Test Method for Protective Clothing Material Resistance to Puncture."