

# Cardinal Health™ Sterilization Wrap

**Your trusted source for ease and assurance**



# Everything you need, wrapped in a single solution

**Every day, you're tasked to do more with less.** To improve productivity with less time, budget and staff. And do it all while you strive to keep your patients safe and improve the quality of care.

**That's why we created a new choice in sterilization wrap:** a full selection of products that are durable, reliable and cost effective.

## All the selection you want

Whatever your clinicians need, now you have a single product line that can take care of it all. Choose the wrap that best meets your goals of maximizing efficiency, safety and quality of care.


### Two Color:

The blue and green layers help make it fast and easy to differentiate, identify and inspect trays

### Single Layer:

Supports sequential wrapping methodologies

- Made from SMS polypropylene, non-woven material for strength and sterility maintenance.
- Designed with a simple, "cross-stitch" bond pattern for clear, visual feedback of breaches that may occur during handling, storage and/or sterilization.
- Available in six basis weights and a wide range of sizes to meet the needs of today's Central Sterile Department.



# Cardinal Health™ Sterilization Wraps

## Two Color

- Made from two layers of SMS — one layer blue, the other green — bonded together along opposing sides
- Delivers the protection of double-wrapping in a single sheet (simultaneous wrapping) — resulting in time and productivity efficiencies when wrapping and opening<sup>1,2</sup>
- Select the color you want — blue or green, to help differentiate and identify instrument trays by type or use

## Single Layer

- Made from one layer of blue SMS
- Supports simultaneous double or sequential wrapping

## Two Color

	Catalog no.	Size/Dimensions	Case qty.		Catalog no.	Case qty.
CH100	CH1G0012	12 x 12 in. (30 x 30cm)	480	CH200	CH110012	1000
	CH1G0015	15 x 15 in. (38 x 38cm)	480		CH110015	1000
	CH1G0018	18 x 18 in. (45 x 45cm)	480		CH110018	1000
	CH1G0020	20 x 20 in. (50 x 50cm)	480		CH110020	1000
	CH1G0024	24 x 24 in. (60 x 60cm)	240		CH110024	500
	CH1G0030	30 x 30 in. (76 x 76cm)	144		CH110030	300
	CH1G0036	36 x 36 in. (91 x 91cm)	144		CH110036	300
	CH1G0040	40 x 40 in. (101 x 101cm)	120		CH110040	250
	CH1G0045	45 x 45 in. (114 x 114cm)	120		CH110045	250
	CH1G0048	48 x 48 in. (121 x 121cm)	120		CH110048	250
	CH1G0054	54 x 54 in. (137 x 137cm)	96		CH110054	100
	n/a	54 x 72 in. (137 x 182cm)	n/a		CH110072	100
CH200	CH2G0012	12 x 12 in. (30 x 30cm)	480	CH210012	1000	
	CH2G0018	18 x 18 in. (45 x 45cm)	480	CH210018	1000	
	CH2G0024	24 x 24 in. (60 x 60cm)	240	CH210024	500	
	CH2G0030	30 x 30 in. (76 x 76cm)	144	CH210030	300	
	CH2G0036	36 x 36 in. (91 x 91cm)	144	CH210036	300	
	CH2G0040	40 x 40 in. (101 x 101cm)	120	CH210040	250	
	CH2G0048	48 x 48 in. (121 x 121cm)	48	CH210048	100	
	CH2G0054	54 x 54 in. (137 x 137cm)	48	CH210054	100	
CH300	CH2G0072	54 x 72 in. (137 x 182cm)	48	CH210072	100	
	CH3G0024	24 x 24 in. (60 x 60cm)	240	CH310024	500	
	CH3G0030	30 x 30 in. (76 x 76cm)	144	CH310030	300	
	CH3G0036	36 x 36 in. (91 x 91cm)	144	CH310036	300	
	CH3G0040	40 x 40 in. (101 x 101cm)	96	CH310040	250	
	CH3G0045	45 x 45 in. (114 x 114cm)	48	CH310045	100	
	CH3G0048	48 x 48 in. (121 x 121cm)	96	CH310048	100	
	CH3G0054	54 x 54 in. (137 x 137cm)	48	CH310054	100	
CH400	CH3G0072	54 x 72 in. (137 x 182cm)	48	CH310072	100	
	n/a	18 x 18 in. (45 x 45cm)	n/a	CH410018	600	
	CH4G0024	24 x 24 in. (60 x 60cm)	120	CH410024	250	
	n/a	30 x 30 in. (76 x 76cm)	n/a	CH410030	200	
	CH4G0036	36 x 36 in. (91 x 91cm)	72	CH410036	150	
	CH4G0040	40 x 40 in. (101 x 101cm)	72	CH410040	150	
	CH4G0045	45 x 45 in. (114 x 114cm)	48	CH410045	100	
	CH4G0048	48 x 48 in. (121 x 121cm)	48	CH410048	50	
CH500	CH4G0054	54 x 54 in. (137 x 137cm)	48	CH410054	50	
	CH4G0072	54 x 72 in. (137 x 182cm)	24	CH410072	50	
	CH5G0018	18 x 18 in. (45 x 45cm)	288	CH510018	600	
	CH5G0024	24 x 24 in. (60 x 60cm)	120	CH510024	250	
	CH5G0030	30 x 30 in. (76 x 76cm)	120	CH510030	250	
	CH5G0036	36 x 36 in. (91 x 91cm)	72	CH510036	150	
	CH5G0045	45 x 45 in. (114 x 114cm)	48	CH510045	100	
	CH5G0048	48 x 48 in. (121 x 121cm)	48	CH510048	50	
	CH5G0054	54 x 54 in. (137 x 137cm)	48	CH510054	50	
	CH5G0060	60 x 60 in. (152 x 152cm)	24	CH510060	50	
CH600	CH5G0072	54 x 72 in. (137 x 182cm)	24	CH510072	50	
	CH5G0090	54 x 90 in. (137 x 228cm)	24	CH510090	50	
	CH6G0036	36 x 36 in. (91 x 91cm)	72	CH610036	144	
	CH6G0040	40 x 40 in. (101 x 101cm)	48	CH610040	96	
	CH6G0045	45 x 45 in. (114 x 114cm)	48	CH610045	96	
	CH6G0048	48 x 48 in. (121 x 121cm)	48	CH610048	48	
CH6G0054	54 x 54 in. (137 x 137cm)	48	CH610054	48		
CH6G0072	54 x 72 in. (137 x 182cm)	24	CH610072	48		

1. Webster, Joan et al. Barrier properties and cost implications of a single versus a double wrap for storing sterile instrument packs. American Journal of Infection Control, Volume 33, Issue 6, 348 - 352.

2. Rutala WA, Weber DJ. Choosing a sterilization wrap. Infection Control Today 2000;4:64,70.

# Cardinal Health™ Sterilization Wraps (continued)

## Paper

- Made from wetlaid cellulose — blue in color
- Single layer format to support sequential or simultaneous double wrapping techniques
- Available in one basis weight and 16 sizes

## Paper

Catalog no.	Size/Dimensions	Case qty.
4008	8 x 8 in. (20 x 20cm)	1000
4009	9 x 9 in. (23 x 23cm)	1000
4012	12 x 12 in. (30 x 30cm)	1000
4015	15 x 15 in. (38 x 38cm)	500
4018	18 x 18 in. (45 x 45cm)	500
4020	20 x 20 in. (50 x 50cm)	500
4024	24 x 24 in. (60 x 60cm)	500
4030	30 x 30 in. (76 x 76cm)	250
4036	36 x 36 in. (91 x 91cm)	250
4040	40 x 40 in. (101 x 101cm)	250
4045	45 x 45 in. (114 x 114cm)	250
4048	48 x 48 in. (121 x 121cm)	100
4054	54 x 54 in. (137 x 137cm)	100
4060	40 x 60 in. (101 x 152cm)	100
4065	60 x 60 in. (152 x 152cm)	100
4072	54 x 72 in. (137 x 182cm)	100

**We are your  
trusted source for  
everything you need**

to be efficient,  
productive —  
and save money  
at the same time.

Ask your Cardinal Health representative  
for a trial today or call customer service  
at 877.CARDINAL.



[cardinalhealth.com](http://cardinalhealth.com)

**Spunbond/Meltblown/Spunbond (SMS):** Consists of meltblown polypropylene between layers of spunbonded polypropylene. It is strong, fluid-resistant, breathable and provides durable protection.

© 2016 Cardinal Health. All Rights Reserved. CARDINAL HEALTH, the Cardinal Health LOGO and ESSENTIAL TO CARE are trademarks or registered trademarks of Cardinal Health.  
Lit. No. 2PERI16-529512 (05/2016) PRC (03/2016)

Cardinal Health  
Waukegan, Illinois