

Woman and Baby



You are the angels.

We are your wings.



CardinalHealth

Essential to care™



**Products
you need.**

**Products
you trust.**

**Now from
one partner —
Cardinal Health.**

Labor and delivery essentials

Intrauterine pressure catheters (IUPCs) _____	2
Fetal spiral electrodes (FSEs) _____	3
Transducer belts _____	4
Abdominal binders _____	5
C-Section straps _____	6
Ultrasound gels _____	6
Lubricating jelly _____	7
Baby board straps _____	7
Infant beanies _____	8
Infant footprinter pads _____	9
Infant tape measure _____	9
Infant heel warmers _____	10
Cord clamp clippers _____	10



Mom and baby care

Mom

Feminine hygiene products _____	12
Perineal cold packs _____	13
Nursing pads _____	13

Baby

Baby wipes _____	14
Baby oil _____	15
Baby lotion _____	15
Baby powder _____	16
Baby shampoo _____	16
Baby shampoo/body wash _____	17
Brush and comb set _____	17

Neonatal intensive care

Umbilical vessel catheters (UVCs), anchors and trays _____	20-21
Peripherally inserted central catheters (PICCs) and trays _____	22
Replogle suction catheters _____	23
Temperature probe covers _____	24
I.V. limb support boards _____	24
Baby shades – phototherapy eye protectors _____	25
Baby adhesive tape _____	25
CPAP nasal cannula kits _____	25
Index _____	26-27



You are the angels.

We are your wings.

Cardinal Health offers a comprehensive portfolio of clinically-trusted devices and products to support the needs of your patients. We are focused on the healthcare needs of women before, during and after pregnancy. We provide quality products to help care for women and keep them comfortable — when they need it the most.

Additionally, our infant care portfolio provides you with quality newborn medical devices and essentials.

When life is precious and new, how can you be sure you're doing all you can to care for your patients?

Start with quality products that are safe and clinically trusted.
Start with Cardinal Health.



Labor and delivery essentials

Intrauterine pressure catheters



Kendall™ IUPC




Transducer tip




Catheter body



Zero/Rezero button



Cable check



Fluid resistant connector

Kendall™ IUPC

Intrauterine pressure catheters are intended to measure uterine activity during labor. The patented Kendall™ IUPC was developed to provide reliable tracings and convenience helping the clinician deliver the best patient care.

- Small, rounded transducer tip allows for easy insertion
- Flat catheter body minimizes trauma by gently deflecting away from tissue during insertion
- Clear, flexible amnioinfusion line delivers 55 mL/mm at gravity feed for fast infusion of liquids
- Convenient indicator light on Zero/Rezero button provides additional confirmation during the rezero process
- Unique Cable Check feature allows clinicians to identify a cable issue
- Fluid resistant connector minimizes erroneous readings by protecting the connection from fluid splash
- Not made with natural rubber latex, DEHP or PVC
- Sterile


Code	Description	Pkg string
56300	Kendall™ IUPC	10 ea/case
Reusable cables		
56311	For most Philips™ monitors	1 ea/case
56312	For Philips™ Avalon™ monitors	1 ea/case
56321-	For GE® Corometrics™ monitors	1 ea/case
56330	For Spacelabs™ monitors	1 ea/case
56313	For Neoventa™ monitors	1 ea/case

IUPC cables



Product sold in case quantities only.

Fetal spiral electrode



Kendall™ FSE Leg plate



Electrode body



FSE attachment pads

Kendall™ FSE System

The FSE system is comprised of three individual parts:

1. A fetal spiral electrode (FSE)
2. An attachment pad
3. A reusable cable

All components are needed in order to successfully use the Kendall FSE System. The fetal spiral electrode is intended to collect fetal heart rate signal and deliver it to a monitor during the delivery process. The attachment pad is used to secure placement of the electrode by adhering to the patient's thigh. Finally the reusable cable is needed to connect the electrode to a fetal monitoring system.

Features

- Protected lead wires with a safety cap
- Spiral is recessed to help reduce risk of needle sticks when removing from package
- Not made with natural rubber latex or DEHP
- Sterile

Benefits

- One step removal of the drive and guide tube
- Flexible guide tube can be gently shaped for easier insertion
- Simple connection to unique leg plate for secure connections, no orientation required
- Attaches to the patient via a self-adhesive grounding electrode

Code	Description	Fetal monitoring system	Pkg string
FSE			
31479549	Kendall™ FSE 2100		50 ea/case
Attachment pads			
31481396	FSE Attachment Pad (small 1.875 x 1.875 in.) disposable, latex-free		50 ea/case
50000095	FSE Attachment Pad (large 3.75 x 1.875 in.) disposable, latex-free		50 ea/case
Reusable cables			
HP®			
31479754	FCB 100	8040A	1 ea/case
Philips™			
31479846	FCB 300	For most Philips™ monitors – 2 ft.	1 ea/case
50000331	FCB 308	For most Philips™ monitors – 8 ft.	1 ea/case
GE®			
31479812	FCB 400	Monitors with round connector	1 ea/case
50000107	FCB 700	Monitors with rectangular connector	1 ea/case

FSE cables



FCB 100



FCB 300



FCB 308



FCB 400



FCB 700

Product sold in case quantities only.

Customer Service: **800.964.5227**

Transducer belts



Kendall™ Abdominal Belts

Kendall™ Abdominal Belts are used to hold toco and/or ultrasound transducers to a mother's abdomen during labor.

- Made from soft materials intended to provide patient comfort
- Large variety of sizes and packaging configurations available
- Retains shape and holds transducers snugly



Buttonhole belt

Code	Description	Pkg string
Buttonhole style		
40000007-	2.375 x 48 in., knit elastic, solid pink and blue set, not made with natural rubber latex	100 sets/case
31410270B-	2.375 x 48 in., knit elastic, solid pink and blue set with 2 bobbins, not made with natural rubber latex	50 sets/case
31410270	2.375 x 48 in., knit elastic, solid pink and blue set, not made with natural rubber latex	50 sets/case
56101	2.375 x 48 in., knit elastic, white with pink and blue stripes, not made with natural rubber latex	50 sets/case



Tab belt

Tab style		
40000009	1.5 x 36 in., knit elastic, solid pink and blue set, not made with natural rubber latex	100 sets/case
40000008	1.5 x 36 in., knit elastic, solid pink and blue set, not made with natural rubber latex	50 sets/case
40000008B	1.5 x 60 in., knit elastic, solid pink and blue set, not made with natural rubber latex	50 sets/case
56102-	1.5 x 42 in., knit elastic, pink and blue striped set, not made with natural rubber latex	50 sets/case
30955689	2 x 52 in., mesh, white, not made with natural rubber latex	50 sets/case



Mesh belt

Product sold in case quantities only.

Abdominal binders



Binders with 2-way stretch

Kendall™ Abdominal Binders with 2-way stretch

Kendall™ Abdominal Binders are used to hold toco and/or ultrasound transducers to a mother's abdomen during labor.

- Vertical and horizontal stretch is designed to provide patient comfort
- Multilateral compression intended to provide secure placement of the transducer
- Breathable material to help wick away moisture
- Finished bands at top and bottom prevent rolling-up of the binder without fraying at the ends
- Color-coded stitching in the waistband quickly identifies sizing
- Individually packaged product
- Not made with natural rubber latex

Code	Description	Color	Width	Length	Pkg string
31312P	MD Pre-cut binder	Blue	12 in.	14 in.	50 ea/case
31314P	LG Pre-cut binder	Black	14 in.	14 in.	50 ea/case
31316P	XL Pre-cut binder	Green	16 in.	14 in.	50 ea/case



Economy binder rolls

Kendall™ Economy Abdominal Binder Rolls

Kendall™ Economy Abdominal Rolls are used to hold toco and/or ultrasound transducers to a mother's abdomen during labor.

- Must be cut to length
- Not made with natural rubber latex

Code	Description	Color	Width	Length	Pkg string
31363362	SM Roll binder	Beige	10 in.	33 ft.	1 roll/case
31364360	MD Roll binder	Beige	12 in.	33 ft.	1 roll/case
31364378	LG Roll binder	Beige	14 in.	33 ft.	1 roll/case

C-section straps



Kendall™ C-Section Straps

Used to secure a patient's legs during a c-section procedure, Kendall™ C-Section Straps are soft, comfortable and not made with natural rubber latex.

Code	Description	Pkg string
30951696	3 x 84 in. C-Section Strap	25 ea/case

Ultrasound gel




Cardinal Health™ Ultrasound Gel

Cardinal Health™ Ultrasound Gel is a salt-free, formaldehyde-free and latex-free formula.

- Acoustically correct
- Dye-free
- Greaseless

Code	Description	Pkg string
USG20ST	Ultrasound gel, 20g, Sterile	48 ea/case
USG20-NS	Ultrasound gel, 20g, Non-sterile	100 ea/bx; 8 bx/case
USG-250BT	Ultrasound gel, 250mL	12 ea/bx; 8 bx/case
USG-5L	Ultrasound gel, 5L	4 ea/case

Lubricating jelly



Cardinal Health™ Lubricating Jelly

Cardinal Health™ Lubricating Jelly helps ease discomfort during a variety of patient examinations and may be used for lubrication of rubber products, such as catheters and rectal thermometers.

- Sterile
- Water soluble
- Greasless
- Unscented
- CHG-free

Code	Description	Pkg string
LJ33121	Lube jelly, 2 oz., flip cap	12 ea/bx; 8 bx/cs
LJ33169	Lube jelly, 4 oz., flip cap	12 ea/bx; 6 bx/cs
LJ33107G	Lube jelly, 3 gm, foil pack 144	144 ea/bx; 6 bx/cs
LJ33183G	Lube jelly, 5 gm, foil pack 144	144 ea/bx; 6 bx/cs
LJ33107	Lube jelly, 3 gm, foil pack	30 ea/bx; 48 bx/cs

Baby board straps



Argyle™ Baby Board Straps

Argyle™ Baby Board Straps gently restrain infants during circumcisions and other routine procedures.

- Not made with natural rubber latex
- Foam
- Nonsterile

Code	Description	Pkg string
31309480	Four 10 in. foam straps per set	4 ea/set; 25 sets/cs
31368858	Four 10 in. straps per set (made of esther/tricot)	4 ea/set; 25 sets/cs

Infant beanies



Argyle™ Infant Beanies

Made from a knitted cotton-polyester blend, our colorful Argyle™ Infant Beanies help to reduce heat loss in newborns and infants. The beanies measure approximately 4 x 5 in.

Available in the following options:

- Single-ply
- Double-ply



Code	Description	Pkg string
Unisex pink and blue striped		
52000010	Argyle™ Double-ply	50 ea/carton; 9 cartons/cs; 450 ea/cs
30953130	Argyle™ Single-ply	50 ea/carton; 9 cartons/cs; 450 ea/cs

Infant footprinter pads



Infant footprinter pads

Our disposable footprinter ink pads will help you create a priceless keepsake for your patients.

Code	Description	Pkg string
IFP-PAD	Infant footprinter pad	25 ea/bx; 4 bx/cs

Infant tape measure



Infant tape measure

Our easy-to-read, disposable, paper tape measure is for single-patient use.

Code	Description	Pkg string
30942B	Medical tape measure infant, 24 in., disposable	500 ea/bg; 20 bg/cs; 10,000 ea/cs

Infant heel warmers



Cardinal Health™ Infant Heel Warmer

Cardinal Health™ Infant Heel Warmers help increase an infant's blood flow so caregivers can get a quality heel stick. Just activate and it's ready to use. They require no cooling off period and have no sharp edges to harm caregiver or patient.

- Consistent and effective heat therapy
- Available in squeeze or click activation
- Heel warmer quickly reaches the proper temperature and maintains this temperature for several minutes
- Not made with natural rubber latex or DEHP

Code	Description	Pkg string
Cardinal Health™ infant heel warmer		
11460-010T	Infant heel warmer squeeze activation	25 ea/bx; 4 bx/cs
11470-010T	Infant heel warmer disc activation	25 ea/bx; 4 bx/cs

Cord clamp clippers



Argyle™ Umbilical Cord Clamp Clippers

Argyle™ Umbilical Cord Clamp Clipper is a disposable, individually packaged, non-sterile device used to remove an umbilical clamp.

- Compatible with most standard hinged clamps
- Not made with natural rubber latex
- Single patient use


Code	Description	Pkg string
52000008	Cord clamp clipper	60 ea/ct; 8 ct/cs



Mom and baby care

Feminine hygiene products

Made to cover the spectrum of absorbency needs



Absorbency



Pads

Code	Description	Individually Wrapped	Absorbency	Width	Length	Pkg string
FH-92144	OB underpad	No	Heavy	9.25 in.	22 in.	24 ea/pk; 10 pk/cs
FH-PADPCH3	Contoured, peach pad	No	Heavy	5.51 in.	13.77 in.	30 ea/pk; 8 pk/cs
FH-PADORG	Ultra thin, overnight maxi pad with wings	Yes	Heavy	6.26 in.	12.2 in.	14 ea/pk; 12 pk/cs
2022A	Curity™ maternity pad	Yes	Heavy	3.89 in.	12.0 in.	14 ea/bg; 12 bg/cs; 168/cs
1580A	Curity™ winged peri-pad	Yes	Moderate	5.86 in.	11.02 in.	16 ea/pk; 12 pk/cs
1380A	Curity™ regular peri-pad	Yes	Moderate	3.6 in.	10.11 in.	24 ea/bg; 12 bg/cs; 288/cs

Perineal cold packs



Cardinal Health™ Perineal Cold Packs

Cardinal Health™ Perineal Cold Packs help provide consistent cold therapy for convenient and effective pain relief. The all-in-one design is an alternative to traditional postpartum temperature therapy and can be applied by the patient.

- Comfortable, absorbent material contains a cold pack that is instantly activated by folding and shaking
- Premium features a liquid-impermeable backing for added comfort and security
- Not made with DEHP

Code	Description	Pkg string
11500-010	Standard 12 1/2 x 4 in.	24 ea/cs
11447-010	Premium 13 1/2 x 6 in.	24 ea/cs

Nursing pads



Curity™ Nursing Pads

Curity™ Nursing Pads are soft and comfortable for nursing care.

- Appropriate for use with regular or nursing bras
- Disposable
- Soft, non-irritating edges and wet-proof backing
- Promotes patient comfort

Code	Description	Pkg string
2630-	Curity™ Nursing Pad, 4.25 x 3.75 in.	12 ea/bx; 24 bx/cs; 288 ea/cs

Baby wipes



Premium wipes

Our thick, plusher wipes come out of the package one at a time, allowing you to stay focused on your newborn — not the product.

All wipes are free of dyes, alcohol, lanolin, parabens, sulfates, 1,4-Dioxane and Q-15 and are not made with phthalates.

- Hypoallergenic and fragrance-free

Code	Description	Pkg string
2BWPU-16	Wipe baby sensitive fragrance-free	16 ea/pk, 48 pk/cs
2BWPU-42	Wipe baby sensitive fragrance-free	42 ea/pk, 24 pk/cs



Standard wipes

Lighter weight fabric, but with a gentle formulation that cleanses naturally. Choose from two formulations: fragrance-free and scented. Both are free of dyes, alcohol, PBA, phthalates, sulfates, 1,4-Dioxane and Q-15.

- **Fragrance-free:** Hypoallergenic with unique fold for easy, 1 wipe at a time dispensing
- **Scented:** Fresh scent with unique fold for easy, 1 wipe at a time dispensing

Code	Description	Pkg string
2BWSU-12	Wipe baby fragrance-free, travel pack	12 ea/pk, 48 pk/cs
2BWSU-40	Wipe baby fragrance-free, flip top	40 ea/pk, 24 pk/cs
2BWSU-80	Wipe baby fragrance-free, tub	80 ea/pk, 12 pk/cs
2BWSS-40	Wipe baby scented, flip top	40 ea/pk, 24 pk/cs

Baby oil



Baby oil

Moisturize baby's skin with 100% mineral oil that creates a barrier on the skin. Formulated without any dyes, fragrances, alcohol, lanolin, parabens, phthalates, sulfates, 1,4-Dioxane or Q-15.

- Free of fragrance and dyes
- Great for baby massages and caregiver hands

Ingredient list: 100% mineral oil

Available in 1 oz. clear, easy-to-use flip-top cap bottle

Code	Description	Pkg string
2BO01	1 oz.	144 ea/cs

Baby lotion



Baby lotion

Our baby lotion gently moisturizes and is hypoallergenic. And of course, it only includes the ingredients needed — such as glycerin and oils to moisturize.

- Free of dyes, fragrance, alcohol, lanolin, parabens, sulfates, 1,4-Dioxane and Q-15
- Not made with phthalates
- Gentle moisturizer
- Hypoallergenic

Available in 1 oz. and 2 oz. clear, easy-to-use flip-top cap bottles

Code	Description	Pkg string
2BL01	1 oz.	144 ea/cs
2BL02	2 oz.	144 ea/cs

Baby powder



Baby powder

Keep baby's skin cool and comfortable, without fragrance or talc. Our baby powder absorbs excess wetness and helps eliminate friction.

- No talc
- No fragrance
- Helps eliminate friction
- Absorbs excess wetness
- Keeps skin cool and comfortable

Ingredient list: cornstarch, sodium bicarbonate, hydrated silica

Available in 1.5 oz. rounded square bottle with dial

Code	Description	Pkg string
2BP15	1.5 oz.	96 ea/cs

Baby shampoo



Baby shampoo

Naturally-derived ingredients are the heart of every baby product. And they're the natural way to protect while you clean. In short, our shampoos are everything you need — and nothing you don't need.

- Free of dyes, fragrance, alcohol, lanolin, parabens, phthalates, sulfates, 1,4-Dioxane and Q-15
- Tear-free formulation
- Naturally low foaming
- Hypoallergenic
- Gentle cleansing

Code	Description	Pkg string
2BS01	1 oz.	144 ea/cs
2BS02	2 oz.	144 ea/cs

Baby shampoo/body wash



Shampoo/body wash

Add the gentle cleansing action of our shampoo and body wash and you have a head-to-toe solution. Like our baby shampoo, this product is also free of unnecessary and harsh chemicals.

- Free of dyes, fragrance, alcohol, lanolin, parabens, phthalates, sulfates, 1,4-Dioxane and Q-15
- Tear-free formulation
- Naturally low foaming
- Hypoallergenic
- Gentle cleansing

Code	Description	Pkg string
2BA01	1 oz.	144 ea/cs
2BA02	2 oz.	144 ea/cs

Brush and comb set



Brush and comb set

Two-piece set, in a gender-neutral ivory. The brush features soft bristles and is packaged with a wide tooth comb.

- Individually wrapped


Code	Description	Pkg string
BBRCMSET	Gender-neutral ivory	24 ea/bx; 12 bx/cs






Neonatal intensive care


Umbilical vessel catheters (UVCs)



Single-lumen with side clamp



Dual-lumen



Triple-lumen

Argyle™ Umbilical Vessel Catheters (UVCs)

This catheter is inserted into the umbilicus and used to administer fluids for hydration, parenteral nutrition, and other common intravenous medications during the first 7 days after birth. UVCs benefit the patient by reducing the need for painful venipunctures and providing double or triple the access of a single lumen catheter. This allows simultaneous therapies from one device.

Features

- Clear, lipid resistant luer lock hubs
- Color coded I.V. extensions on triple and dual lumen catheters
- Sterile
- Not made with DEHP
- Not made with natural rubber latex

Benefits

- Polyurethane material for increased maximum strength and flow rates
- Encapsulated X-ray opaque lines
- Numerical depth markings every 1 cm from 5-25 cm
- Effective length of 25 cm

Code	French size	Primary/Secondary/ Tertiary Lumen size	Length	Pkg string
Polyurethane single-lumen UVC with luer lock hubs				
8888160325	2.5	26ga	11.5 in. (29.21 cm)	10 ea/cs
8888160333	3.5	21ga	14 in. (35.56 cm)	10 ea/cs
8888160341	5.0	18ga	13.6 in. (34.544 cm)	10 ea/cs
Polyurethane dual-lumen UVC with luer lock hubs				
8888160531	3.5	20ga / 23ga	15 in. (38.1 cm)	5 ea/cs
8888160556	5	18ga / 21ga	15 in. (38.1 cm)	5 ea/cs
Polyurethane triple-lumen UVC with luer lock hubs				
8888160648	5	18ga / 21ga / 21ga	15 in. (38.1 cm)	5 ea/cs
8888160663	8	15ga / 18ga / 18ga	15 in. (38.1 cm)	5 ea/cs

Umbilical catheter anchor



Argyle™ Umbilical Catheter Anchor

Provides a fast, convenient, strain-relief bridge used to hold UVCs in place.

- With its two adhesive arms, the umbilical catheter anchor can hold two separate lines securely in place

Code	Description	Pkg string
MI30000	Umbilical catheter anchor	50 ea/cs

Umbilical vessel catheter insertion trays



Argyle™ UVC Insertion Tray

Argyle™ Umbilical Vessel Catheter Insertion Trays

Contents

PREP TRAY (top section)

- 1 - Umbilical vessel catheter (optional) in pre-packaged sterile catheter system
- 1 - Pack povidone iodine swab sticks
- 1 - Fenestrated drape
- 2 - Towels
- 6 - Gauze pads 2 x 2 in.

PROCEDURE TRAY (bottom section)

- 1 - CSR wrap
- 2 - Filter aspiration needles
- 1 - 3-way stopcock
- 1 - Stainless steel scissors
- 1 - Stainless steel needle holder, 5 in.
- 1 - Stainless steel eye dressing forceps, half-curved
- 1 - Stainless steel eye dressing forceps, fully curved
- 2 - Stainless steel hemostats, curved
- 1 - Stainless steel hemostat, straight
- 1 - Blade scalpel #11
- 1 - Stainless steel Adson tissue forceps
- 1 - Stainless steel vessel probe dilator
- 1 - Monoject™ 1cc syringe
- 2 - Monoject™ 6cc syringes
- 1 - Monoject™ 3cc syringe
- 1 - 3.0 Silk suture with curved needle
- 1 - Umbilical tape 15 in.
- 1 - Syringe tip cap
- 6 - Gauze pads 4 x 4 in.
- 1 - Monoject™ 20ga blunt needle**
- 1 - Monoject™ 18ga blunt needle**
- 1 - Safety scalpel #11***
- 2 - Monoject™ 20ga x 1.5 safety needle***

**Only included in trays without catheter - 8888160408 and 160409

***Only included in tray - 43201

Code	Description	Pkg string
8888160119	UVC tray with 3.5 and 5 fr. polyurethane UVC with luer lock hub	5 ea/cs
8888160424	UVC tray with 3.5 fr. polyurethane UVC	5 ea/cs
8888160432	UVC tray with 5 fr. polyurethane UVC	5 ea/cs
8888160838	UVC tray with 3.5 fr. dual lumen polyurethane UVC	5 ea/cs
8888160853	UVC tray with 5 fr. dual lumen polyurethane UVC	5 ea/cs
8888160408	UVC tray without catheter and with fine tip forceps	5 ea/cs
160409	UVC tray without catheter and with smooth tip forceps	5 ea/cs
43201	UVC tray without catheter and with safety components	5 ea/cs

Neonatal-pediatric peripherally inserted central catheters (PICCs)



1.9 Fr dual lumen PICC

Argyle™ Neonatal-Pediatric Peripherally Inserted Central Catheters (PICCs) and Trays

This catheter may be used to administer fluids for hydration and parenteral nutrition, as well as other common intravenous medications, in case where long term I.V. administration is necessary.

Features

- Not made with natural rubber latex or DEHP
- Easy-to-read graduation marks
- Stabilizing wings
- Clear pigtails for visualization of fluid
- 30 cm catheter length
- 0.6 mm outer diameter

Benefits


- Polyurethane material for increased maximum strength and flow rates
- Radiopaque material for verification of proper placement



Procedure tray

The dual lumen PICC has all of the above features and benefits, PLUS

- Independent lumens throughout the entire catheter to help prevent mixing of incompatible infusates
- Staggered pigtails to differentiate primary and secondary lumens and for ease of securement and maintenance



1.9 Fr single lumen PICC

Code	Description	Pkg string
Single lumen products		
43303	1.9 Fr single lumen PICC, not made with natural rubber latex	10 ea/cs
43309	1.9 Fr single lumen PICC procedure tray	5 ea/cs
Dual lumen products		
43304	1.9 Fr dual lumen PICC, not made with natural rubber latex	10 ea/cs
43311	1.9 Fr dual lumen PICC procedure tray	5 ea/cs
Procedure tray only		
43308	PICC procedure tray, no catheter	5 ea/cs
Introducer		
43320	Peel away cannula introducer	5 ea/cs

PICC tray contents:

- 1 - Poly lined drape
- 1 - Measuring tape

Prep tray

- 1 - Measuring tape
- 1 - Tourniquet
- 1 - Fenestrated drape
- 1 - Towel
- 3 - 2 in. x 2 in. Gauze
- 6 - 4 in. x 4 in. Gauze
- 1 - Pack of povidone iodine swab sticks
- 1 - Pack of alcohol swab sticks


Procedure tray

- 1 - Wound closure strips
- 1 - 12cc syringe
- 1 - MONOJECT™ MAGELLAN™ 18ga X 1.5 in. safety needle
- 1 - Scissors
- 2 - Iris forceps, half-curved
- 1 - Transparent dressing



Peel away cannula introducer

Replogle suction catheters



Argyle™ Replogle Suction Catheters

The Argyle™ Replogle Suction Catheter is used to treat Esophageal Atresia in neonatal patients.

- Not made with DEHP
- Clear dual-lumen tube
- Sterile, three smooth eyes
- X-ray opaque Sentinel Line™ tubing with Sentinel Eye™ catheter
- Not made with natural rubber latex
- Made with PVC material

Code	Description	French Size	Length	Pkg String
8888256545	6FR Replogle suction catheter	6	24 in. (61 cm)	10/case
8888256529	8FR Replogle suction catheter	8	24 in. (61 cm)	10/case
8888256503	10FR Replogle suction catheter	10	24 in. (61 cm)	10/case

Temperature probe covers



Argyle™ Temperature Probe Covers

Our disposable temperature probe covers are designed with the delicate skin of your smallest patients in mind.

- Gentle hydrogel adhesive securely anchors temperature probe to infant's skin
- May be removed and repositioned
- Not made with natural rubber latex

Code	Description	Pkg string
Hydrogel		
31424784	Duck and bear shaped (small)	200/case
31424792	Duck and bear shaped (large)	200/case
MI30002	Rectangular (.75 x 1.25 in.)	800/case

I.V. limb support boards



Argyle™ I.V. Limb Support Boards

The Argyle™ I.V. Limb Support Boards are specially designed, pre-shaped boards that promote immobilization of the patient's limb while the I.V. is in place.

- Repositionable--can be bent and re-bent while the patient's I.V. is in place
- Not made with natural rubber latex
- Optional VELCRO™ strap provides help with securing to the patient's limb
- Fabric and foam cover to promote patient comfort

Code	Description	Pkg string
MI00880	1 x 3 in. (preemie)	50/case
MI00775	2 x 4 in. (small)	50/case
MI00850	2 x 5 in. (large)	50/case

Baby shades – phototherapy eye protectors



Argyle™ Baby Phototherapy Eye Protectors

The Argyle™ Baby Phototherapy Eye Protectors are used to protect the eyes during the treatment of infantile jaundice.

- Adheres to infant's temples with our biocompatible hydrogel adhesive or optional strap
- Not made with natural rubber latex

Code	Description	Pkg string
MI00575-	1.125 x 4.5 in. (2.8 x 10.8 cm) (small) for infants less than 1500 grams	50 ea/carton; 4 cartons/cs; 200 ea/cs
MI00625	1.5 x 5.5 in. (3.4 x 13.9 cm) (large) for infants more than 1500 grams	50 ea/carton; 4 cartons/cs; 200 ea/cs

Baby adhesive tape



Argyle™ Baby Adhesive Tape

The Argyle™ Baby Adhesive Tape is designed with our special clear hydrogel formula to gently yet firmly secure catheters, I.V. lines, lead wires and other medical devices to the fragile, sensitive skin of infants.

- Clear adhesive hydrogel formula gently yet firmly secures catheters, I.V. lines and lead wires to infant's skin
- Tapes may be repositioned and reapplied without irritating or damaging skin
- Not made with natural rubber latex

Code	Description	Pkg string
MI00677	.75 x 2 in. tape strips	50 strips/pack; 10 pk/case

CPAP nasal cannula kits



Argyle™ CPAP Nasal Cannula Kits

The Argyle™ CPAP Nasal Cannula is a sterile, single-use disposable oxygen delivery system with rigid plastic connector and soft-molded silicone cannula with two short-angled prongs.

- Silicone material
- Hook and loop straps
- Carefully angled nares
- 15mm connectors
- Not made with natural rubber latex

Code	Description	Size	Pkg string
8888163055	For babies under 1000 grams	X-Small	10/case
8888163030	For babies 1000-1500 grams	Small	10/case
8888163014	For babies over 1500 grams	Large	10/case

Index

CODE	PRODUCT	PAGE
------	---------	------

Abdominal binders

31312P	MD Pre-cut binder	5
31314P	LG Pre-cut binder	5
31316P	XL Pre-cut binder	5
31363362	Economy SM roll binder	5
31364360	Economy MD roll binder	5
31364378	Economy LG roll binder	5

Baby adhesive tape

MI00677	.75 x 2 in. tape strips	25
---------	-------------------------	----

Baby board straps

31309480	Four 10 in. foam straps per set	7
31368858	Four 10 in. straps per set	7

Baby lotion

2BL01	1 oz.	15
2BL02	2 oz.	15

Baby oil

2BO01	1 oz.	15
-------	-------	----

Baby powder

2BP15	1.5 oz.	16
-------	---------	----

Baby shades – phototherapy eye protectors

MI00575-	Small	25
MI00625	Large	25

Baby shampoo

2BS01	1 oz.	16
2BS02	2 oz.	16

Baby shampoo/body wash

2BA01	1 oz.	17
2BA02	2 oz.	17

Baby wipes

2BWPU-16	Premium sensitive	14
2BWPU-42	Premium sensitive	14
2BWS-40	Standard scented	14
2BWSU-12	Standard fragrance-free, travel pack	14
2BWSU-40	Standard fragrance-free	14
2BWSU-80	Standard fragrance-free, tub	14

Brush and comb set

BBRCMSET	Gender-neutral ivory	17
----------	----------------------	----

Cord clamp clippers

52000008	Cord clamp clipper	10
----------	--------------------	----

CODE	PRODUCT	PAGE
------	---------	------

CPAP nasal cannula kits

8888163014	Large	25
8888163030	Small	25
8888163055	X-Small	25

C-Section straps

30951696	C-section strap	6
----------	-----------------	---

Feminine hygiene products

1380A	Curity™ Regular Peri-Pad singles	12
1580A	Curity™ Winged Peri-Pad singles	12
2022A	Curity™ Maternity Pad, unscented	12
FH-92144	OB underpad	12
FH-PADORG	Unscented, ultra thin, overnight maxi pad	12
FH-PADPCH3	Unscented, contoured, peach pad	12

Fetal spiral electrodes

31479549	Kendall™ FSE 2100	3
31479754	Philips™ FCB 100 8040A	3
31479812	GE FCB 400 Monitors with round connector	3
31479846	FCB 300 for most Philips™ monitors–2 ft	3
31481396	FCB 2000 (small 1.875 in.)	3
50000095	FCB 2001 (large 3.75 in. x 1.875 in.)	3
50000107	GE FCB 700 Monitors with rect. connector	3
50000331	FCB 308 for most Philips™ monitors–8 ft.	3

Infant beanies

30953130	Argyle™ Single-Ply	8
52000010	Argyle™ Double-Ply	8

Infant footprinter pads

IFP-PAD	Infant footprinter pad	9
---------	------------------------	---

Infant heel warmers

11460-010T	Infant heel warmer squeeze activation	10
11470-010T	Infant heel warmer disc activation	10

Infant tape measure

30942B	24 in.	9
--------	--------	---

Intrauterine pressure catheters

56300	Kendall™ IUPC	2
56311	Reusable – for most Philips™ monitors	2
56312	Reusable – for Philips™ Avalon™ monitors	2
56313	Reusable – for Neoventa™ monitors	2
56321-	Reusable – for GE Corometrics™ monitors	2
56330	Reusable – for Spacelabs™ monitors	2

I.V. limb support boards

MI00775	2 x 4 in. (small)	24
MI00850	2 x 5 in. (large)	24
MI00880	1 x 3 in. (preemie)	24

CODE	PRODUCT	PAGE
------	---------	------

Lubricating jelly

LJ33121	Lube jelly, 2 oz., flip cap _____	7
LJ33107G	Lube jelly, 3 gm, foil pack 144 _____	7
LJ33169	Lube jelly, 4 oz., flip cap _____	7
LJ33183G	Lube jelly, 5 gm, foil pack 144 _____	7
LJ33107	Lube jelly, 3 gm, foil pack _____	7

Neonatal-pediatric peripherally inserted central catheters (PICCs)

43303	1.9 Fr single lumen PICC _____	22
43304	1.9 Fr dual lumen PICC _____	22
43308	PICC procedure tray, no catheter _____	22
43309	1.9 Fr single lumen PICC procedure tray _____	22
43311	1.9 Fr dual lumen PICC procedure tray _____	22
43320	Peel away cannula introducer _____	22

Nursing pads

2630-	Curity™ Nursing Pad, 12 _____	13
-------	-------------------------------	----

Perineal cold packs

11447-010	Premium 13 1/2 x 6 in. _____	13
11500-010	Standard 12 1/2 x 4 in. _____	13

Replogle suction catheters

8888256503	Replogle suction catheters, 10 _____	23
8888256529	Replogle suction catheters, 8 _____	23
8888256545	Replogle suction catheters, 6 _____	23

Temperature probe covers

31424784	Duck and bear shaped (small) _____	24
31424792	Duck and bear shaped (large) _____	24
MI30002	Rectangular (.75 x 1.25 in.) _____	24

Transducer belts

56101	Buttonhole style _____	4
30952090	Buttonhole style _____	4
31410270B-	Buttonhole style _____	4
31410270	Buttonhole style _____	4
40000007-	Buttonhole style _____	4
56102-	Tab style _____	4
30955689	Tab style _____	4
40000008B	Tab style _____	4
40000008	Tab style _____	4
40000009	Tab style _____	4

Ultrasound gels

USG-5L	Ultrasound gel, 5L _____	6
USG20-NS	Ultrasound gel, 20g NS _____	6
USG20ST	Ultrasound gel, 20g STL _____	6
USG-250BT	Ultrasound gel, 250mL _____	6

CODE	PRODUCT	PAGE
------	---------	------

Umbilical catheter anchor

MI30000	Umbilical catheter anchor _____	20
---------	---------------------------------	----

Umbilical vessel catheter insertion trays

43201	UVC tray w/o cath. and w/ safety components _____	21
160409	UVC tray w/o cath. and w/ smooth tip forceps _____	21
8888160119	3.5 and 5 fr. polyurethane UVC w/ luer lock hub _____	21
8888160408	UVC tray w/o catheter and w/ fine tip forceps _____	21
8888160424	3.5 fr. polyurethane UVC _____	21
8888160432	5 fr. polyurethane UVC _____	21
8888160838	3.5 fr. dual lumen polyurethane UVC _____	21
8888160853	5 fr. dual lumen polyurethane UVC _____	21

Umbilical vessel catheters (UVCs)

8888160325	Single-lumen UVC 2.5 _____	20
8888160333	Single-lumen UVC 3.5 _____	20
8888160341	Single-lumen UVC 5.0 _____	20
8888160531	Dual-lumen UVC 3.5 _____	20
8888160556	Dual-lumen UVC 5 _____	20
8888160648	Triple-lumen UVC 5 _____	20
8888160663	Triple-lumen UVC 8 _____	20

Just the right products for the right outcomes.

For more information or to order, contact your Cardinal Health sales representative today, call **877.CARDINAL** or visit **cardinalhealth.com**.

© 2020 Cardinal Health. All Rights Reserved. CARDINAL HEALTH, the Cardinal Health LOGO, CURITY, KENDALL, ARGYLE, MONOJECT, MAGELLAN, SENTINEL EYE and SENTINEL LINE are trademarks or registered trademarks of Cardinal Health and may be registered in the US and/or in other countries. All other marks are the property of their respective owners. Patent www.cardinalhealth.com/patents. Lit. No. 2NP20-1191420 (6/2020)



cardinalhealth.com