

## Upholstered Mat Platform

### MODELS: 240-47/240-57/240-68

- Heavy duty, upholstered mat platform in Natural wood tone
- 6-leg design with center support adds stability
- Triple-bolted, corner legs
- Solid hardwood legs and plywood frame
- Frame attached with steel angles
- 2" ultra-firm, foam padding
- Gray edge bumper guard, helps protect top from tears
- Heavyweight, knit-backed vinyl upholstery
- Optional 21" total height available – see options
- 700 lbs. load capacity under normal use



California residents please be advised of the following, pursuant to Proposition 65:

**WARNING:** This product can expose you to chemicals, including diisononyl phthalate (DINP), which is known to the State of California to cause cancer. For more information go to [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).

2

#### MODEL 240-47 SPECIFICATIONS

Top 7' L x 4' W x 18" H

#### MODEL 240-57 SPECIFICATIONS

Top 7' L x 5' W x 18" H

#### MODEL 240-68 SPECIFICATIONS

Top 8' L x 6' W x 18" H

### Options



041 Maximum 21" High

### Standard Upholstery



Slate Blue

Wedgewood

Royal Blue

Desert Tan

Burgundy

Country Mist

Warm Gray

Gunmetal

Black