



Avalo Anesthesia Cart/AM10MC-EG-K-DR321

Brand: Capsa Healthcare

Item#: CPAM10MC-EG-K-DR321

Material: High-density polyethylene

Dimensions: W31" x H43" x D24"

- * Six drawers: (3) 3", (2) 6", (1) 9"
- * Keyless Entry lock
- * Keyless Entry with Auto Relock available
- * Proprietary "Pyramid" shape side panel provides lower center of gravity to prevent tipping
- * Precision-welded steel frame with high-impact cart panels
- * Seamless drawer with full extension for maximum access
- * Secure "drawer detent" closure
- * Slot Lock Advantage - All drawers within 1" of closed will be secured when cart is locked

Alternative heights and drawer configurations available

Additional colors available

Anesthesia Cart Accessory Package (ACP-AM10MC-EG-1) Includes:

- * Anesthesia storage tower
- * 4 and 5 front tilt bins
- * Anesthesia tray
- * 3" drawer divider kit
- * Sharps case
- * Glove dispenser
- * Accessory mounting rack
- * Two handles

Additional accessories available

--