

Life Begins at the End of the Ramp!!

Don't be stopped by the edges of ramps, or beach mats!

DeBug Beach Wheelchair gives you the freedom to enjoy the beach, or any outdoor terrain the way you deserve.

Without Restrictions! The DeBug is designed by the Demings who know first hand the obstacles wheelchair users face on a daily basis. They have used their personal experience and incorpated it nto the design and quality of their products.

LONG LASTING #316L STAINLESS STEEL FRAME / ADJUST-ABLE TELESCOPIC LEG REST / FITS THROUGH A STANDARD ADA DOORWAY / LOW MAINTENANCE DELRIN BUSHINGS / POLYURETHANE PNEUMATIC WHEELS ROLL ACROSS SAND, SNOW AND GRAVEL / DUAL STITCHED SLING WITH TENERA 130 THREAD / PHIFERTEX PLUS SLING MATERIAL RESISTS STRETCHING, MILDEW AND BACTERIA, COLD CRACK AND ULTRAVIOLET DAMAGE AND IS MACHINE WASHABLE.



SEAT SIZES	ENGLISH (INCH)	METRIC (CM)
Widths	14, 18, 22	36, 46, 56
Average depth	17	43
Seat-to-floor height	17 - 21	43 - 53
Backrest height	18	46
Telescope footrest knee to heel	15 - 22	38 - 56
Elevating leg rest knee to heel	15 - 22	38 - 56

WEIGHTS Capacity 250 - 400 lbs 114 - 182 kg Chair weight 40 - 65 lbs 18 - 30 kg

debug







ARMREST STYLES

4 - 182 kg Removable swing away armrest

18 - 30 kg Adjustable height or padded optional

"Contoured" or specialized optional

Interchangeable desk option available