

# Specifications

**De-Bug frame** is constructed using #316L Stainless Steel tubing which is a higher grade than standard “#304 marine grade” stainless.

- \* **Weight capacity** 350 pounds standard, weight capacities up to 600 pounds can be accommodated for a nominal fee.
- \* **Casters** - rear wheels rotate 360 degrees . Wheel hubs are made from poly-nylon material to eliminate salt water corrosion.
- \* The rear suspension also articulates +/- ~20 degrees, allowing all four wheels to remain in contact with the ground when going across uneven terrain, making it much more stable than standard outdoor wheelchairs.
- \* **Swing away arm rests** are also easily removable for Lateral transfers.
- \* **Footrest**-Has Linear position telescopic length **adjustment** and four position angle adjustment. The De-Bug also has a footrest for proper foot support.
- \* **Wheel Lock** - Adjustable tension wheel brake can be actuated independantly on each side, uses all #316L stainless steel hardware.
- \* **Chair weight** - Forty pounds for a standard adult chair
- \* **Chair width** is ~33.5 inches to access a 36 inch doorway and is ~27 inches wide with the optional interchangeable pool wheels.
- \* **Chair design** is patented in the U.S. #380,991

## 30 DAY Money Back Guarantee

### Ten year warranty

Components & Frame #316L series Stainless Steel “the higher the number the higher quality of Stainless”

### Two year warranty

Seating - Material is Phifertex Nylon, anti-absorbant, U.V. resistant. Seating is removable and machine washable.

### Two year warranty

**Tires** - 3 ply poly urathene, low pressure, to specifically enable easy operation in grass, gravel, sand, snow, & uneven terrain



*Deming Designs Inc.  
1090 Cobblestone Dr.  
Pensacola, FL 32514  
850-478-5765  
850-476-3361 Fax  
[www.beachwheelchair.com](http://www.beachwheelchair.com)*

