

Deluxe Aluminum Shower Bench without Back

HCPCS: E0245



FIGURE A



FIGURE B



FIGURE C



FIGURE D



FIGURE E

- Tool-free assembly of seat and legs (Figure A)
- Drainage holes in bench reduce slipping
- Aluminum frame is lightweight, durable and corrosion proof
- Legs are height adjustable in 1" increments, and angled with suction-style tips to provide additional stability
- Made using recycled materials

PRODUCTS	Description
RTL12203KDR	Bath Seat w/o Back, KD, Retail,1 EA
SPECIFICATIONS	Description
Warranty	Limited Lifetime
Seat Dimensions	19.25" (W) x 11.5" (D)
Weight Capacity	400 lbs.
Carton Shipping Weight	7 lbs
Weight	5 lbs 5 lbs 4.5 lbs