

## Deluxe Automatic Blood Pressure Monitor, Upper Arm

HCPCS: A4670-1,A4670



FIGURE A



FIGURE B

- WHO Classification
- Irregular Heart Beat Detection
- One Touch Operation
- X-Large Easy-to-Read Screen
- Date & Time Display
- Pulse Rate Indicator
- Digital Error Messages
- Travel Bag
- Auto-Off Function
- Memory stores 120 readings (2 x 60)
- Average Last 3 Results

PRODUCTS	Description
BP3400	BP Monitor, Auto, Upper Arm, 1/EA
SPECIFICATIONS	Description
Auto Power Off	3 Minutes After Last Reading
Carton Size	18.9" x 17.7" x 17.7"
Cuff Size	8.5"Ø14"
LCD Size	2.45" x 1.81"
Measurement Method	Oscillometric Method
Measurement Range	Pressure: 0mmHg~300 mmHg Pressure: +/- 3 mmHg
Memory Capacity	2 x 60 Reading with Date & Time
Operation Environment	50°F to 104°F, <85% RH
Power Supply	4x AA Batteries or AC Adapter
Unit Dimensions	5.83" x 3.94" x 2.2"
Warranty	2 Years
ACCESSORIES	Description
BP2001	Cuff,Std,BP3400,BP3600,BPHL,8.5-14, EA
BP2000	Cuff,PlusSize,12.5"x17.75",BP3400,3600, EA
BP2610	A/C Adapter, BP3400,BP3600,1/ea, EA