



The 553KL is perfect for newborns, infants and toddlers. The scale operates on batteries or AC power so you can use it anywhere. Practices benefit from one-button operation that makes it easy to use while ensuring that the baby remains safely on the scale.

Advantages

- EMR Connectivity**
Scale has the ability to interface with a computer, EMR software or other electronic device via various Health o meter[®] Professional Connectivity Solutions
- Value-Added Features**
Operating functions (LB/KG Conversion, KG Lock Out, Zero Out/Tare, Hold/Release, Auto Off) included at no additional cost
- Motion Detection**
Scale has the ability to weigh accurately despite movement of the subject on the scale
- Measuring Tape Included**
Easy-to-read measuring tape allows user to measure both weight and height simultaneously



Carrying Case
Item # 553CASE



Specifications

- Capacity: 44 lb / 20 kg
- Resolution: 0.5 oz / 10 g
- Tray Size (w x d x h):
24 7/8" x 13" x 1 3/8"
632 mm x 330 mm x 35 mm
- Product Footprint (w x d x h):
24 7/8" x 16 1/2" x 3 7/8"
632 mm x 419 mm x 98 mm
- Product Weight: 10 lb / 5 kg
- Measuring Tape: 0" - 24" / 0 cm - 58 cm
- Display: 1" / 25 mm LCD
- Power Source: 6 AA batteries included, optional power adapter (order # ADPT30)
- Connectivity: USB, optional Pelstar wireless technology
- Carrying Case: Optional carrying case (order # 553CASE)
- Functions: Zero, Tare, Hold / Release, Auto Off
- 2 Year Limited Warranty

Model	Description
553KL	Digital pediatric tray scale, pounds & kilograms
553KG	Digital pediatric tray scale, kilograms only
553KL-HR	553KL with attachable mechanical height rod
553KL-EHR	553KL with attachable digital height rod
SS-553KL	ScaleSurance 2 year extended warranty for 553KL & 553KG*
553HR	Attachable mechanical height rod for 553KL
553EHR	Attachable digital height rod for 553KL
2210CART	Rolling cart for 553KL
553CASE	Carrying case for 553KL
ADPT30	Optional 100-240V adapter
PLUGSET30	Plug adapters for use outside North America
PELSTARONE-P	Wireless thermal printer