

Digital Physician's Scale

KEY FEATURES

For today's busy physicians' offices, Doran's DS5100 Digital Physician's Scale provides a faster and easier way to weigh and measure patients accurately. With our *Touchless Start* feature, you can acquire a patient's weight in as little as two seconds without any buttons to press. Simply have the patient step onto the large platform and an accurate weight up to 500 pounds will be held on the display. Once the patient steps off the device, the DS5100 is immediately ready for another weighing.

The DS5100's sturdy aluminum column houses the easy-to-use height measurement rod. Simply extend the measurement arm and pull upwards for adult heights or push down to measure a child. To calculate BMI, use the scroll buttons to enter the height, then press the BMI button and the patient's *Body Mass Index* is instantly displayed in the right-hand LCD. When not in use, the BMI display conveniently functions as a clock in 24 or 12-hour format.

The DS5100 is EMR Ready and data can be transmitted via the included USB port.

SPECIFICATIONS

Warranty:	Three Years
Power Input:	Four C Batteries (Included)
Battery Life:	20,000 Weighments
Maximum Capacity:	500 x 0.1 lb (200 x 0.05kg)
Accuracy:	+/- 0.1 lb (0.05kg)
Dimensions:	15" W x 13" D x 2" H
Height Rod Adjustment:	Minimum 2' - 7" (81 cm) Maximum 6' - 10" (210 cm)
Display:	Dual 1" High LCD
Displayed Units:	lb, kg
Controls:	Zero, LB, KG, Height Input/BMI, FT/IN, CM, Scroll Up/Down and Clock
Clock:	12 or 24 Hour Display
BMI Calculator:	Input data in FT/IN or CM
Construction:	ABS and Metal Construction
Communications:	EMR Ready with USB Standard
Option	USB Communication Cable

