

PreserveTech Aquachair Bathing System with Bidet

HCPCS:



FIGURE A



FIGURE B



FIGURE C




FIGURE D

Enjoy thorough, independent bathing with the PreserveTech Aquachair Bathing System with Bidet. The perfect solution for individuals with unsteadiness, standing difficulties or problems reaching, the premium shower seat features a convenient in-seat bidet plus handheld shower sprayer and an open seat to help the user clean hard-to-reach areas for optimal hygiene. For additional security and easier transfers, add the optional PreserveTech Aquachair Transfer Accessory with Curtain Guard Protection (Item # RTL12A004-TR).

- Bring Confidence Back to Bathing
- Unique shower chair with bidet and shower spray enables full-body cleansing without standing
- Choose front or rear cleansing modes (front cleansing is gentler than the rear spray)
- Easily access shower spray from armrest holder. Open seat makes it easy to access hard-to-reach areas.
- Padded non-slip seat and backrest offer support and comfort
- Adjust the height with the push of a button to achieve the ideal seated position. Seat height adjustments on this model range from X" to X"
- Legs include a non-skid, non-slip tip for additional security on wet surfaces
- Lightweight and durable-rust-resistant aluminum frame supports up to 350 lbs.
- Anti-microbial treated backrest, seat and tips help prevent stains and odors from buildup of bacteria, mold and mildew.
- Easily connects to the water supply no electrical connection or batteries necessary

PRODUCTS

Description

RTL12A004-WH	AQ Shower Chair & Commode, 1EA, 1CV
--------------	-------------------------------------

SPECIFICATIONS

Description

Warranty	1yr Limited
Product Weight	10.89 Lbs
Material	PP + Aluminium
Weight Capacity	350 lbs
Color	White
Dimensions Assembled	32.28" (W) X 18.11" (L) X 34.65" (H)
Retail Box Dimensions	12.2" (W) X 17.52" (L) X 21.65" (W)
Retail Box Weight	2.43 lbs
Width Between Arms	17.8"

CA Prop 65	yes WARNING: This product can expose you to chemicals, including Styrene, which is known to the State of California to cause cancer. For more information go to www.P65Warnings.ca.gov .
------------	--

ACCESSORIES

Description

RTL12A004-TR	AQ Transfer Bench Accessory, 1EA, 1CV, EA
--------------	---