

CHAD® Evolution Electronic Oxygen Conserver

HCPCS:



FIGURE A

- Delivers a uniform oxygen pulse at increments: 1, 2, 3, 4, 5, 6, 7 Plus CF preset @ 2 LPM
- Minimum 2-year battery life* utilizing 2 AA batteries
- Delivers a 5:1 savings ratio from 14-40 breaths per minute
- Senses a breath/delivers oxygen in the first half of the inspiratory cycle
- Operating pressure range: 500-3000 PSI
- Weighs only 14.9 ounces with batteries
- Single push button selects the equivalent LPM dosage settings
- Auto shut-off feature conserves battery life
- Fits easily in virtually any carry bag

PRODUCTS

Description

OM-900	Evolution Conserver,1/ea
--------	--------------------------

SPECIFICATIONS

Description

Warranty	Limited 2 Year
----------	----------------

ACCESSORIES

Description

SEAL100	Seal Washer,Conservers/Regulators,10/pkg, PKG
OM-900-BD	BatteryDoors&Screws,OM-900/OM-900M,5/pkg, EA