

Forearm Platforms for Nimbo

HCPCS: E1399-1



FIGURE A

- Constructed of lightweight aluminum
- Height and depth adjustable
- Handgrips extend forward or back, supinate, pronate and are angle adjustable
- Contoured armrest with a flexible edge
- Includes hook and loop fastener straps to prevent arms from slipping
- Choice of mounting positions - inside or outside of the handlebar for lateral adjustment Spacers are included for additional lateral range

PRODUCTS

Description

KA 1035 FPL	Nimbo Forearm Platfrm/MtgBrckt,Lg,1pr/bx
KA 1035 FPS	Nimbo Forearm Platfrm/MtgBrckt,Sm,1pr/bx

SPECIFICATIONS

Description

Armrest Length	5" 8"
----------------	---------

PARTS

Description

CE 0041	Ratchet/Bolt,CE1035&CE1045, 1/ea, EA
KA 1035 FPPL	Mounting Plate, Large, KA 1035 FPS / FPL, EA
KA 1035 FPHDW	Hardware for KA 1035 FPS / FPL, 1/ea, EA
KA 1035 KNOB	Knob, KA 1035 FPS / FPL, 1/ea, EA
CE 1035 EC	EndCap,vert post,CE1035,1/ea, EA
KA 1035 FPS	Nimbo Forearm Platfrm/MtgBrckt,Sm,1pr/bx, BX
KA 1035 FPL	Nimbo Forearm Platfrm/MtgBrckt,Lg,1pr/bx, BX