

Overbed Table, Double Top



HCPCS:



FIGURE A

- Features two table surfaces - -a larger one that tilts in either direction (Figure A) and a smaller one that stays securely in place
- Easily assembled
- Walnut, wood grain low-pressure laminate top
- Tops can be raised or lowered in infinite settings between 30" - 46"
- Swivel casters allow for easy maneuverability
- Chrome-plated steel "H" base provides security and stability
- Tabletops lock securely when height adjustment handle is released
- Tabletop level can be raised with the light upward pressure

PRODUCTS**Description**

13068BV	OVER BED TABLE, DOUBLE TOPS, BV, 1/ea
---------	---------------------------------------

SPECIFICATIONS**Description**

Base Dimensions	25" (W) x 28" (D)
Carton Shipping Weight	24.5 lbs
Height Adjustment	30" - 46"
Tabletop Surface Dimensions	22.8" (W) x 15" (D) (Larger) 15" (W) x 15" (D) (Smaller)
Warranty	Limited Lifetime
Weight	22 lbs
Weight Capacity	55 lbs

PARTS**Description**

13007CL	Caster for 13007 w/Lock,1/ea, EA
13007C	Caster,Overbed table.1/ea, EA
13068BV-AAR	AngleAdjustmentRodW/Bushing&Knob,13068BV, EA
13068BV-AKST	AngleAdjustmntKnob,SmallTableTop,13068BV, EA
13068BV-BRKTLLT	BracketForLargeTableTopW/Hrdwr,13068BV, EA
13068BV-LTT	Large Table Top, 13068BV, 1/ea, EA
13068BV-MAST	Mast, 13068BV, 1/ea, EA
13068BV-HDWR	Hardware Pack, 13068BV, 1/ea, EA
13068BV-STT	Small Table Top, 13068BV, 1/ea, EA
13068BV-BRKT-STT	Bracket For Small Table Top,13068BV,1/ea, EA