

Aluminum Rollator, 7.5" Casters

HCPCS: E0143+E0156



FIGURE A



FIGURE B



FIGURE C



FIGURE D



FIGURE E



FIGURE F



FIGURE G




FIGURE H

- Comes with a comfortable, durable padded seat
- Zippered pouch under seat for added privacy and security of personal items
- Removable, hinged, padded backrest can be folded up or down
- 7.5" casters with soft grip tires ideal for indoor and outdoor use and make the ride smooth and stable
- Ergonomic handles are easy-to-grip and are adjustable

PRODUCTS	Description
R728BL	Rollator,7.5" 4Whl,Alum,Blue,EA
R728RD	Rollator,7.5" 4Whl,Alum,Red,EA

SPECIFICATIONS	Description
Casters	7.5"
Dimensions	24W x 28D
Handle Height	33" - 38"
Inside Hand Grip Width	17"
Material	Aluminum
Number of Wheels	4
Seat Dimensions	14" (W) x 12" (D)
Seat To Floor Height	23"
Warranty	Limited Lifetime
Weight	14 lbs
Weight Capacity	300 lbs

PARTS	Description
9505W1026120	Tote,10261,726,728,R800,R900, 1/ea, EA
10208KN	Knob,Ht.Adj.Rollator,w/Attch.1/ea, EA
9501S726BK	Back Rest with Foam,Black,726,728, 1/ea, EA
9500001232	Brake Cable/Adjust,Lower,1/ea, EA
750CLIP	Wheel LockClip,750N,1/ea, EA
9500004101	Cable,Fastner,1/ea, EA
9501S72800-01	Front Wheel for 728/R728,1/ea, EA
9502M72801	FrontFork/CasterAssembly,728,2M Ser#,1ea, EA
950340059901-2M	Padded Seat Brackets, 2M Ser.#, 1/pr, PR
9501S72603	Backrest Pivot Brackets,R726,10257,1/ pr, PR
9501S72616	Spring,R726,1/ea, EA

PARTS

Description

9500075011	Rear Axle Bolt,750,780,1/ea, EA
9501S7261L	Wheel Lock, Left, R726,1/ea, EA
9501S7261R	Wheel Lock, Right, R726,1/ea, EA
9501S72800	Rear Wheel for 728/R728,1/ea, EA
STDS1002-1	Axle for 8" caster, 1/ea, EA
9501S72602	Hand Grip for 726,1pr, PR
9501S72609	Folding Hinge Assembly for 726/728,1/ea, EA
9500002211	Axle Bearing,Rollator,1ea, EA
10208C	CableClip for RollatorBrake,1/ea, EA