

Airial Spotz® the Dog Mask with Meter Dose Inhaler Chamber

HCPCS:




FIGURE A

- Momentarily holds aerosol medication in chamber until patient is ready to inhale
- Accommodates most CFC and HFA meter dose inhalers
- Helps to eliminate waste of medication
- Minimizes need for hand-breath coordination
- Includes Spotz(R) pediatric mask for use with chamber
- Maximizes drug delivery
- Easy-to-use and clean
- Anti-static properties and Latex Free

PRODUCTS

Description

MQ8100

Airial Chamber with Spotz Mask 1/ea