

## Pulmo-Aide® Compressor Nebulizer System

HCPCS: E0570,A7003



- The Pulmo-Aide Compressor Nebulizer System is the only approved compressor for use with TOBI (Tobramycin)
- One of the longest lasting compressors on the market with an average product life of 10 years
- Unmatched quality and durability Great for nursing homes and long-term users
- Metal frame motor construction
- Large built-in compartment for storing nebulizers and tubing
- Built-in nebulizer holder is a convenient place to store your nebulizer if your treatment is interrupted

PRODUCTS	Description
5650D	Pulmo-Aide w/ Disp. Nebulizer
SPECIFICATIONS	Description
Warranty	Limited 5 Year
Dimensions	6.3" (W) x 10.5" (H) x 10.1" (D)
Free Air Flow	12 lpm Typical
Weight	7.1 lbs
ACCESSORIES	Description
5650D-601	FILTER PKG - DFT (25-PK), PKG
5650D-602	FILTER PKG - MM, 5/PK, PKG
5650D-636	COMPRESSOR TEST FIXTURE, EA
5650D-614	POWER SWITCH (6/PK), EA
5650D-604	REAR COVER ASSY, EA
5650D-603	FRONT COVER ASSY, EA
5650D-605	DOOR ASSY, EA
5650D-609	PUSH-IN FEET (4-PK), EA
5650D-611	TUBE CONNECTOR W/NUT-6-PK, PKG
5650D-612	FILTER KIT(CAP&FILTER) (6/PK), PKG
DL1050	PEDIATRIC CHARACTER MASK (50-PK), CS
3655D-621	VixOne Disp. Nebulizer, CS 50, CS