



Deluxe Folding Travel Walker with 5" Wheels

HCPCS: E0143









FIGURE A

FIGURE B

FIGURE C

FIGURE D





FIGURE E

FIGURE F

- Easily folds tool-free to ½ the size of a traditional walker (Figure A B) and can be stored in a "Ready Set Go" Pouch (Figure C)
- Legs are attached with strong elastic cords so legs can be folded (Figure D)
- Comes with 5" front wheels and rear glide skis that are ideal for most surfaces
- Two side-storage pouches with large pockets for personal items are held in place by hook and loop fastener tabs to both side frames

PRODUCTS	Description
RTL10263KDR	Walker w/ 5" Whls,Travel,Red,1 EA
SPECIFICATIONS	Description
Carton Shipping Weight	10 lbs.
Height	29.5 - 37"
Overall Width	26.5"
Warranty	Limited Lifetime
Weight	8 lbs.
Weight Capacity	300 lbs
Width	17" inside hand grip 22" inside back legs