

Shower Chair with Back and Removable Padded Arms

HCPCS:



FIGURE A

- Removable padded armrests provide comfort and added security
- Blow-molded bench and back provide comfort and strength
- Drainage holes in bench reduce slipping
- Aluminum frame is lightweight, durable and corrosion proof
- Leg height adjusts in 1" increments
- Made using recycled materials

PRODUCTS	Description
12445KD-1	KD Bath Bench w/Back & Padded Arm,1/ea

SPECIFICATIONS	Description
Carton Shipping Weight	10 lbs. 21 lbs.
Outside Legs	17.75 (W) x 19.25 (D)
Seat Depth	14.5", 16.25"
Seat Height	16" - 21"
Seat Width	16"
Warranty	Limited Lifetime
Weight	9.75 lbs. 19 lbs.
Weight Capacity	300 lbs
Width Between Arms	18.5"

ACCESSORIES	Description
12002T	Tips for 12002,4/pkg., BX