

## Saratoga Spring Bathing System













FIGURE A

FIGURE B

FIGURE C

FIGURE D



FIGURE E

- Thermostatic mixing valve Provides an anti-scald protection that dispenses 15 gallons per minute for fast tub fill
- Built-in disinfection system The most time-efficient solution for disinfection of the baths external surfaces and minimizes the risk of cross-infection
- The built-in disinfection system cabinet is conveniently located on the back end of the tub for easy access
- The lockable cabinet provides a secure place for cleaning supplies and disinfectants
- · Handheld shower wand Hair and upper body washing can be used at the same time the tub is filling
- · Water temperature display Monitors tub-fill and shower wand temperatures
- 3/4" supply lines Constructed of stainless steel mesh clad
- Three plumbing configurations to choose from side, front, or island configuration
- The island configuration allows 360° unobstructed access for the caregiver and with no visible hook-ups or cords, the Saratoga Springs Bathing System will give the bathing area the look of a luxury spa
- The CleanDrive Air Spa System Variable speeds and pulse modes designed to enhance the bathing experience for the resident and caregiver

Description
Reclining Tub w Thermo Valve, Disinfect
Reclining Tub w Thermo Valve, Disinfect
Description
24" (W) x 21" (H)
450 lbs or Resident 730 lbs or Resident and Water
78"
34" (W) x 73" (L) x 51" (H)