

Travelite Transport Chair



HCPCS:



FIGURE A



FIGURE B




FIGURE C



FIGURE D

- The Travelite Chair in a Bag is a lightweight yet sturdy aluminum transport chair in gunmetal gray
- It is easily transportable and comes with a compact carry bag with wheels making it ideal for traveling (Figure A - B)
- Folds for storage and transportation (Figure C)
- It is ideal for use indoors and outdoors
- 18" square nylon seat comes with lift-up straps
- Comes with wheel locks 8" solid tires
- Comes with padded, flip-up desk arms and height adjustable footrests

| PRODUCTS | Description |
|--------------------------------|---|
| TC005GY | Travelite Chair in a Bag,1/ea |
| SPECIFICATIONS | Description |
| Back Height | 17.5" |
| Carton Shipping Weight | 23.2 lbs |
| Overall Dimensions | 22.5" (W) (open) 11.5" (W) (closed) x 37.5" (H) |
| Overall Length (with riggings) | 33" |
| Seat Dimensions | 18" (W) x 18" (D) |
| Seat To Floor Height | 19" |
| Warranty | Limited Lifetime |
| Weight | 20.5 lbs |
| Weight Capacity | 250 lbs |
| PARTS | Description |
| TC00503 | Upholstery,TC005,1ea, EA |
| TC00504L | ArmRestAssemblyLeftCompleteTC005,1ea, EA |
| TC00509 | Foot Rest Connector,TC005,1ea, EA |
| TC005012 | Foot Rest Assembly,TC005,1pair, PR |
| TC005027 | Lifting Strap for TC005GY, 1/ea, EA |
| TC00507 | Arm Rest Release Knob,TC005,1ea, EA |
| TC005011 | CrossFrameBlocksSetof4TC005,1package, PKG |
| STDS1060 | Wheelchair Fork Bearing,1/ea, EA |
| TC005020 | FrontForkTC005,1ea, EA |
| STDS4S3602 | Caster bearing, 1/pr, PR |
| TC005021F | FrontWheelWithBearingsTC005,1ea, EA |
| STDS1011 | Wheel Bearing, 8mm x 22mm 1/ea, EA |
| STDS3J2422 | Frame Cap,Cruiser3,1/ea, EA |

PARTS

Description

TC005024

Carry Bag,TC005,1/ea, EA