

Patient Interfaces

DreamEasy™ Nasal Mask

The DreamEasy™ Nasal CPAP Mask provides a superior seal with its removable Comfort Cushion™ that conforms around the nose and facial contours for a comfortable night's rest. Equipped with a 360° rotational elbow and standard 22 mm tubing connector allowing for freedom of movement when sleeping. Comfortable forehead silicone pad diffuses any pressure. Quick-release clips allow for easy removal of the mask without readjusting the headgear. Dishwasher safe for easy cleaning. Replacement Comfort Cushion™ seals available. New Starter Kit includes all three sizes and a sizing guide for new users.




Product Specifications

- Extremely lightweight and durable
- Flexible design allows for freedom of movement while maintaining a secure fit
- Headgear with Velcro® fasteners provide maximum versatility and movement
- 360° rotational elbow with standard 22 mm tubing connector
- The Comfort Cushion™ made with medical liquid silicone is completely odorless
- Comfort Cushion™ can be removed for easy cleaning
- Ergonomic design is airtight
- Vent hole provides effective CO2 washout
- Latex-free
- Mask is dishwasher safe for easy cleaning
- 90 day limited warranty
- HCPCS Code - A7034* (mask)
- HCPCS Code - A7032* (replacement seal)
- HCPCS Code - A7035* (headgear)



Replacement Seal



Sizing Guide
(CPM-DENSK Only)



Roscoe Medical Offers a 30-Day Money-Back Guarantee on all Roscoe® CPAP Masks, effective from the date of purchase.

If your patient is not 100% completely satisfied and it is within 30 days from your Roscoe purchase date, we will refund the purchase price to you, no questions asked.

Item #	Description	Qty
CPM-DENS	DreamEasy™ Nasal Mask with Headgear, small	each
CPM-DENM	DreamEasy™ Nasal Mask with Headgear, medium	each
CPM-DENL	DreamEasy™ Nasal Mask with Headgear, large	each
CPM-DENSK	DreamEasy™ Nasal Mask Starter Kit with Headgear, All Sizes	each
90091	DreamEasy™ Nasal Mask Replacement Seal, small	each
90093	DreamEasy™ Nasal Mask Replacement Seal, medium	each
90095	DreamEasy™ Nasal Mask Replacement Seal, large	each
90395	Replacement Headgear Clips	2/pack
90377	Replacement Headgear	each

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.



To Order:
 (800) 3-ROSCOE
 fax: (440) 572-4261
 roscoemedical.com