

## Fingertip Pulse Oximeter



HCPCS:



FIGURE A



FIGURE B

- Large easy-to-read multi-view LED display screen
- Displays SpO and pulse rate
- Low power consumption, Auto-off function
- Approximately 30 hours of normal operation
- Accommodates a wide range of finger sizes from pediatric to adult
- Includes lanyard and two AAA batteries.
- Perfect for sports enthusiasts to measure oxygen saturation and pulse rate

### PRODUCTS

	Description
MQ3000	Fingertip Pulse Oximeter 1/ea

### SPECIFICATIONS

	Description
Battery Duration	20 to 30 hours
Battery Type	2 AAA alkaline or rechargeable
Dimensions	2.2(L) x 1.2(W) x 1.3(H)
Warranty	Limited 2 Year
Weight	1.8 oz