



EarthLite®
WORLD'S #1 BRAND IN MASSAGE™

ELLORA™

THE BEST VALUE IN SPA

The Ellora™ line is by far the most popular lift table in the market, it offers a lifetime of use at an affordable price. This heavy duty ADA Compliant table features a solid steel frame with a reliable motor that holds up to 600 lbs. Earthlite's Ellora™ provides you with the perfect balance between reliability and comfort, with deluxe cushioning and upholstery it provides years of trouble-free use.



HEALING HUMANITY THROUGH TOUCH™

ELLORA COLLECTION

Features

- Heavy duty steel frame
- Reliable actuator
- Pro-Plush™ Deluxe 3", 3-layer cushioning system
- Durable and comfortable Comfort-Flex™ upholstery in 6 standard colors: Black, Mystic Blue, Latte, Sterling, Marie's Beige, Vanilla Crème
- Hands free foot control
- Wheels for ease of movement
- ADA compliant
- Lifetime Limited Warranty
- 3 years Warranty on foam and upholstery
- 2 years Warranty on electric parts

Options

- Manual tilt top
- Non-electric pneumatic tilt or salon top
- Electric tilt or salon top with hand control
- Extra foot pedal
- Additional Comfort-Flex™ color selections
- Ultraleather® or Promessa® upholstery
- Flex arms, neck roll and salon pillow
- Flex-Rest™ or Caress™ platform
- Conforma Lux™ Adaptive Chest Comfort System
- 4" Memory foam system
- ETL option

Specifications

- Widths: 28", 30", 32"
- Length: 73"
- Height range: 18" – 37"
- Recommended lift capacity: 600 lbs.
- CE



Flat




Tilt



Salon

Pro-Plush™ Deluxe 3-Layer, 3" Cushioning



Non-Electric Pneumatic Adjustment



Manual Tilt Adjustment

