



## BLADES

# SAFESHIELD STERILE DISPOSABLE SCALPELS

Designed with an ergonomically rounded SUREgrip™ plastic handle that fits the hand comfortably, and outfitted with a stainless steel blade, Safeshield™ Safety Scalpels feature a locking retractable protective shield that helps prevent injuries before, during and after use.

Feather, a blade manufacturer of global renown, utilizes high-precision technology and carefully selected materials. Each Feather scalpel is made in Japan and meets strict ISO standards, resulting in exceptionally high quality and consistent strength and sharpness.

### Features

- **Single Use/Disposable:** Helps prevent patient cross-contamination
- **Individually packaged:** Non-woven Tyvek cloth package complete with clear window for product selection verification and pull tabs for easy opening
- **Clear retractable safety shield:** Can be operated with one hand
- **Auditory confirmation:** Transparent blade shield makes audible click when locked in position with blade covered or uncovered
- **Available in 7 sizes:** Including most common blade sizes and select specialty blade sizes
- **Color-coded by size:** Release tab on protective shield and coordinating packaging assist with blade identification
- **Ruler printed on handle:** Includes cm and mm markings
- Similar weight as conventional disposable scalpel
- **Sterile:** Gamma-sterilized and ready for use immediately after opening
- Individually packaged 10 per dispensing box



Item No.	Stainless Steel Blades	Blade Size	Color Code	UOM
2980#10		10	Red	10/box dispenser
2980#11		11	Pink	10/box dispenser
2980#14		14	Red	10/box dispenser
2980#15		15	Purple	10/box dispenser
2980#20		20	Yellow	10/box dispenser
2980#21		21	Green	10/box dispenser
2980#22		22	Black	10/box dispenser

### Reasons to Make the Switch to Safety Scalpels

- Industry requirements and professional organizations demand them
- Protect employees from on-the-job injuries
- Create a safe place to work; have a vested interest in avoiding the health, staffing, and legal costs of sharp injuries

**NOTE:** Dispose of scalpel in an approved sharps container in accordance with your facility's sharps disposal procedures



GF Health Products, Inc. is an  
intertek ISO 13485:2016 / MDSAP Certified Company

© 2021, GF Health Products, Inc. All Rights Reserved.  
GF Health Products, Inc. is the US distributor for FEATHER  
FEATHER is a registered trademark of FEATHER Safety Razor Co., Ltd



+1 770.368.4700

Information contained herein is subject to change without notice.  
The most current and complete product information can be found on our website.  
[www.grahamfield.com](http://www.grahamfield.com)

Color variations may occur due to monitor or printer settings.

**grahamfield**

GF-1500081-RevD22