

pediatric seating

first class™ school chair



mobile chair shown with optional support kit

- use in classroom, clinic or home
- for children with mild to moderate involvement
- multiple adjustments and accessories meet the needs of each child
- seat is height and depth adjustable with a 15° anterior or posterior tilt
- armrests are height and width adjustable
- pelvic adjustable seat belt helps maintain proper hip position
- optional **support kit** contains H-harness for anterior chest support, adjustable length abductor, individually adjustable height & width lateral supports

chairs	small		large	
stationary chair	31-3800	490.00	31-3802	530.00
mobile chair	31-3810	607.50	31-3812	645.00
accessories				
tray	31-3832	155.00	31-3842	122.50
mobility base	31-3830	205.00	31-3840	160.00
extension legs	31-3838	165.00	31-3838	165.00
support kit	31-3834	300.00	31-3844	277.50
footrest	31-3835	130.00	31-3845	120.00
headrest	31-3833	180.00	31-3833	180.00
hip guides	31-3836	110.00	31-3836	110.00
anti-tip supports	31-3837	50.00	31-3847	45.00

chair dimensions all lengths in inches

	seat width	seat depth	seat height to floor	seat height to footplate	back height	seat to armrest
small	12	9-14	9-13	6-12	15	1.5-9.5
large	16	13-17	14-23	6-22	17	1.5-11.5
	between armrest	between hip guides	between laterals	headrest add'l height	weight capacity	
small	9.5-14.5	8-13	6-13	4.5-8	100 lb	
large	13.5-18.5	12-17	10-17	4.5-8	200 lb	

roll chair

- padded roll seat helps to maintain proper hip abduction
- adjustable angle front & rear legs
- includes lap belt
- small (to 50" tall) - seat 8" dia x 18"D; seat height 8"-16"; seat back 10"W x 14"H
- large (to 68" tall) - seat 10" dia x 22"D; seat height 10"-18"; seat back 12"W x 16"H



roll chair	small		large	
roll chair	31-1120	575.00	31-1121	600.00
accessories				
adjustable clear tray	31-1122	225.00	31-1123	225.00
adjustable lateral supports	31-1124	115.00	31-1124	115.00

sling back chair

- adjustable seat, sling-style adjustable depth back and adjustable footboard
- includes lap belt
- seat 12"W; seat height from floor 14"; depth 6"-12"
- footboard adjusts 5"-11" from seat
- optional tray is wood grain with raised rim; 24"W x 24"D; 7" x 10" cut out



sling back chair		
31-1130	sling back chair	300.00
accessory		
31-1131	optional tray	125.00

multi-use chairs



- stabilizing skis prevent tipping
- adjustable back, footrest and seat depth
- includes lap belt
- small - seat 11½"W x 14"D; depth 6"-12"; seat height from floor 14"; armrest height to floor 18"; footrest adjusts 5"-11" from seat
- large - seat 14"W x 16"D; depth 9"-15"; seat height from floor 17"; armrest height to floor 24"; footrest adjusts 3"-13" from seat



multi-use chair	small		large	
multi-use chair	31-1140	425.00	31-1141	550.00
accessories				
A adjustable clear tray	31-1122	225.00	31-1123	225.00
B adjustable lateral supports	31-1124	115.00	31-1124	115.00
C knee abductor	31-1145	80.00	31-1145	80.00
D adjustable pelvic support	31-1146	140.00	31-1147	150.00
E adjustable foot bracket	31-1148	115.00	31-1149	125.00