

# pediatric positioning

**SKILLBUILDERS®**

## Skillbuilders® strap wedge



- two 5" wide positioning straps provide support when positioning the child
- without the removable straps, the deluxe wedge functions just like the standard Skillbuilders® wedge
- used with optional leg abductor wedge (30-1018 shown on page 219)

	L x W x H		incline	
30-1020	22" x 20" x 4"	small	11°	290.00
30-1021	22" x 20" x 6"	small	16°	305.00
30-1022	22" x 20" x 8"	small	20°	330.00
30-1023	22" x 20" x 10"	small	25°	360.00
30-1024	26" x 24" x 6"	medium	9°	350.00
30-1025	26" x 24" x 8"	medium	13°	370.00
30-1026	26" x 24" x 10"	medium	17°	410.00
30-1027	26" x 24" x 12"	medium	22°	450.00

## Skillbuilders® adjustable wedge system



- arms and hands can move freely while lying on this adjustable wedge
- both the front and rear legs can be adjusted to the proper height with single button locks
- adjustable: front 2-15", rear 2-7"
- removable abductor adjusts 4-8"
- side rail for carrying, securing positioning straps
- wedge folds flat
- vinyl covered foam padded top

31-1050 adjustable wedge (22"W x 36"L) 300.00

## Skillbuilders® bi-height wedge



- use at higher or lower height (side 1 or side 2)
- high molded sides give lateral support
- positioning strap provides support
- accommodates abductor wedge (30-1018 shown on page 219)

	L x W	bi-height wedge				
		side 1	incline	side 2	incline	
30-1240	22" x 18"	4"	11°	6"	16°	310.00
30-1241	22" x 18"	8"	20°	10"	25°	345.00
30-1242	26" x 22"	6"	9°	8"	13°	390.00
30-1243	26" x 22"	10"	17°	12"	22°	440.00

## Versa-Form® positioning pillows



- provides customized, semi-permanent (regular: 1 month, plus: 3 months) positioning
- shape/reshape in seconds by opening air valve and molding to desired shape and then evacuating air with vacuum pump

flat size	Versa-Form®	Versa-Form® Plus	covers
16 x 20"	30-2024 152.50	30-2004 162.50	30-2034 75.00
22 x 12"	30-2020 152.50	30-2000 162.50	30-2030 75.00
22 x 26"	30-2021 162.50	30-2001 172.50	30-2031 77.50
22 x 34"	30-2022 192.50	30-2002 202.50	30-2032 80.00
32 x 40"	30-2023 335.00	30-2003 375.00	30-2033 90.00

## vacuum pump

- single handed operation allows the user to form the pillow

30-2040 pump 122.50

