

Baseline[®] economy SF caliper



with floating measurement tips

- plastic skinfold caliper is the inexpensive choice for skinfold analysis
- spring-loaded arm assures reproducible results
- floating tips assure accurate and reproducible results
- 50mm scale
- complete with booklet and tables

12-1112	economy caliper (each)	20.00
12-1112-25	25 each	450.00



Baseline[®] skinfold calipers

- new 2-sided caliper is easy to use and read for left- or right-handed people; includes calibration block
- lightweight aluminum caliper
- "floating tips" to assure accurate and reproducible results
- large, easy-to-read 70mm scale
- includes protective carrying case
- calibration block verifies at 5-50mm

skinfold calipers


12-1110	skinfold caliper (single sided)	200.00
12-1130	deluxe (dual sided) with calibration block	250.00

accessory

12-1116	calibration block only	25.00
---------	------------------------	-------



single sided



2-sided deluxe skinfold caliper viewable on front and back

double sided

calibration block

hand held BF analyzer



- measures body fat (4% - 45%) and water (37.8% - 66%) with 0.1% resolution
- 8-person memory
- adult and child modes
- lithium batteries included

12-1122	hand-held BF analyzer	30.00
---------	-----------------------	-------

digital FatTrack[®]

- caliper stores 3 measurements, calculates and displays body fat percentage
- Pro model has 50-person memory and carrying case
- includes CR2032 battery




FatTrack™ Pro

12-1124	FatTrack [®] Gold	50.00
12-1129	FatTrack [®] Pro	120.00

Baseline[®] slim guide SF

- for screening large populations
- 60mm scale
- spring-loaded arm
- with instructions and tables



12-1125	slim guide caliper (each)	30.00
12-1125-25	25 each	675.00

BodyLogic[®] analyzer



- body fat measurement in seconds
- no need to pinch skinfolds; just grip the lightweight analyzer with both hands
- stores 9 personal profiles
- five-button keypad

12-1121	BodyLogic	60.00
---------	-----------	-------

Baseline[®] body fat analyzer



- measures/monitors body fat %
- 5 - 50% body fat, 0.1% resolution
- 8-person memory
- 3V battery included
- portable unit, credit card style

12-1141	BF analyzer	10.00
12-1141-25	25 each	225.00

Baseline[®] mini BF analyzer



- measures/monitors body fat %
- 5 - 50% body fat, 0.1% resolution
- 8-person memory
- 3V battery included
- features clock and alarm

12-1140	palm size	10.00
12-1140-25	25 each	225.00

Skyndex[®] skinfold caliper

- automatically calculates body fat %
- digital skinfold caliper available with pre-programmed formulas
- 1 dial for age and gender selection



single formula

12-1100	formula Jackson-Pollock	500.00
12-1101	Durnin	500.00

dual formula

12-1103	Slaughter-Lohmann	500.00
12-1107	Durnin/ Jackson-Pollock	600.00
12-1108	Durnin/ Slaughter-Lohmann	600.00
12-1109	Jackson-Pollock/ Slaughter-Lohmann	600.00