



CanDo® 6" black composite roller

- made of molded foam for increased firmness
- ends don't degrade with use
- purchase by the case and save (in red below)


	round		
	dia x L		
30-2280	6" x 36"	22.50	
30-2280-12	6" x 36"	257.50	
30-2282	6" x 18"	15.75	
30-2281	6" x 12"	11.25	
30-2281-36	6" x 12"	385.00	


	half-round		
	dia x L		
30-2290	6" x 36"	15.75	
30-2290-24	*6" x 36"	360.00	
30-2292	6" x 18"	11.00	
30-2291	6" x 12"	7.90	
30-2291-72	*6" x 12"	540.00	*buy by the case and save



CanDo® 6" EVA foam roller

- rollers are made from extra firm EVA foam
- firm foam maintains soft feel
- attractive mottled design


	round		
	dia x L		
30-2200	6" x 36"	50.00	
30-2201	6" x 12"	25.00	


	half-round		
	dia x L		
30-2210	6" x 36"	35.00	
30-2211	6" x 12"	17.50	



CanDo® 6" blue roller

- blue foam rollers use same industry standard polyethylene foam as white foam rollers
- attractive blue color keeps roller looking new


	round		
	dia x L		
30-2150	6" x 36"	25.00	
30-2151	6" x 12"	12.50	


	half-round		
	dia x L		
30-2152	6" x 36"	17.50	
30-2153	6" x 12"	8.75	




CanDo® Plus antimicrobial foam roller

- extruded polyethylene foam made with anti-microbial additive to protect against microbes and odors
- brown tint differentiates anti-microbial foam from standard foam
- extruded half round has skin on flat side to prevent warping (full-skin)


	round		
	dia x L		
30-2330	6" x 36"	30.00	
30-2331	6" x 12"	15.00	


	half-round		
	dia x L		
30-2332	6" x 36"	21.00	
30-2333	6" x 12"	10.50	



Tufcoat® roller

- TufCoat® rollers have coated exterior that is durable and adds an extra layer of firmness; can be wiped clean
- interior is standard white extruded foam


	round		
	dia x L		
30-2235	4" x 36"	30.00	
30-2236	4" x 12"	15.00	


	half-round		
	dia x L		
30-2245	4" x 36"	21.00	
30-2246	4" x 12"	10.50	



CanDo® full-skin extruded half-round roller

- extruded as a half circle, not sliced in half, providing a full skin
- full skin keeps the half-round roller from warping or bowing and keeps surface cleaner

	half-round		
	dia x L		
30-2340	6" x 36"	21.00	
30-2341	6" x 12"	10.50	



full-skin provides protection and stability on the top and bottom