



FUTURE HEALTH CONCEPTS, INC.

specializing in new and remanufactured medical equipment

Standard Features


- Stainless steel or glass door
- Data logging of Real Time & Temp w/ USB Retrieval
- Keyed temperature lock out
- Push button digital operation
- Large dual digital displays
- Displays actual temperature and set point at all times
- Separate digitally controlled heating chambers
- Rapid warm time (2-6 hours depending on load)
- Uniform heating
- Fully insulated
- Audible and visual over-temp alarm



(USB Retrieval Shown)

Options

- 18" depth
- Mobile base
- Full recess
- Pass through chamber
- Left hand door hinge
- Keyed door lock
- Seismic brackets
- 220/240 V, 50/60Hz
- Custom sizes
- Solution baskets
- Touch screen



SETTING THE INDUSTRY STANDARD