Guide to 3pp[®] Finger Splints Web Pages

Everything you need to know to about 3-Point Products Finger Splints can be found in the following web pages on our website (www.3pointproducts.com), including product details, testimonials, instructional videos, FAQ's, blog posts and where to buy information.

3pp® Buddy Loops®

www.3pointproducts.com/buddy-loop/



Secure, reusable finger protection without adhesives or rough edges.

3pp® Finger Trapper®

www.3pointproducts.com/finger-trapper/



No more fingernail hooks, tape or glue. Use for post-op or dynamic splinting.

Polycentric Hinged Ulnar Deviation Splint™

www.3pointproducts.com/polycentric-hinged/



Reduces ulnar deviation caused by Rheumatoid Arthritis.

3pp[®] Step Up[™] Splint

www.3pointproducts.com/3pp-step-up-splint/



Provides adjustable tension to restore full IP extension.

3pp[®] Step Down[™] Splint

www.3pointproducts.com/3pp-step-down-splint/



Puts the wearer in control to restore IP flexion.

Gel Tubes

www.3pointproducts.com/gel-tubes/



Thin, silicone tubes cushion fingers to relieve pain.

3pp® Final Flexion Wrap®

www.3pointproducts.com/final-flexion-wrap/



Easy-to-apply wrap restores end range IP joint flexion.

Radial Hinged Ulnar Deviation Splint

www.3pointproducts.com/radial-hinged-ulnar -deviation-splint/



Improve alignment of strongly-deviated fingers.

3pp[®] Side Step[™] Splint

www.3pointproducts.com/3pp-side-step-splint/



Reduce lateral deviation at the DIP or PIP joints.

Oval-8® Finger Splints

www.3pointproducts.com/oval-8-finger-splint--copy-1/



Stabilize and align the small joints of the fingers affected by arthritis or injury.



3pp® Finger Splints Sizing & Ordering Information

3pp® Buddy Loops®

www.3pointproducts.com/buddy-loop/

Three sizes. 5" long loops fit most fingers. Sold in packages of 3, 5, 25, 50 or 100.

	3/Pk	5/Pk	25/Pk	50/Pk	100/Pk
1/2" (1.27cm)	RP1005-3	P1005-5	P1005-25	P1005-50	P1005-100
3/4" (1.9cm)	RP1003-3	P1003-5	P1003-25	P1003-50	P1003-100
1" (2.5cm)	RP1007-3	P1007-5	P1007-25	P1007-50	P1007-100

3pp® Finger Trapper®

www.3pointproducts.com/finger-trapper/

One size fits most.

Size	Item Number
Pack of 5	P1006-5
Pack of 25	P1006-25

Polycentric Hinged Ulnar Deviation Splint™

www.3pointproducts.com/polycentric-hinged/

Instructions: Lay hand flat and measure width across the hand from the 2nd to 5th MP (large knuckles).

	Right	Left	Fits
xs	P2002-R1	P2002-L1	2 1/4" to 2 3/4"
Sm	P2002-R2	P2002-L2	2 3/4" to 3"
Med	P2002-R3	P2002-L3	3" to 3 3/8"
Lrg	P2002-R4	P2002-L4	3 3/8" to 3 5/8"

3pp[®] Step Up[™] Splint

www.3pointproducts.com/3pp-step-up-splint/

Instructions: Measure length from web space to fingertip and width at proximal phalanx. Fits right or left hand. If between sizes, choose larger size.

	Right or Left	Web Space to Tip
Sm	P1200-2	1 1/2" - 2" (3.8-5cm)
Med	P1200-3	1 3/4" - 3 1/8" (4.4-8.3cm)
Lrg	P1200-4	3" - 3 1/2" (7.6-8.9cm)

3pp[®] Step Down[™] Splint

www.3pointproducts.com/3pp-step-down-splint/

Instructions: Measure length from web space to fingertip and width at proximal phalanx. Fits right or left hand. If between sizes, choose larger size.

	Right or Left	Web Space to Tip
Sm	P1202-2	1 1/2" - 2" (3.8-5cm)
Med	P1202-3	1 3/4" - 3 1/8" (4.4-8.3cm)
Lrg	P1202-4	3" - 3 1/2" (7.6-8.9cm)

Gel Tubes

www.3pointproducts.com/gel-tubes/

One size fits most.

Size	Item Number
Pack of 5	P1111-5

3pp® Final Flexion Wrap®

www.3pointproducts.com/final-flexion-wrap/

One size fits most.

Size	Item Number
Pack of 5	P1004-5
Pack of 25	P1004-25

Radial Hinged Ulnar Deviation Splint

www.3pointproducts.com/radial-hinged-ulnar-deviation-splint/

Instructions: Lay hand flat and measure width across the hand from the 2nd to 5th MP (large knuckles).

	Right	Left	Fits
Sm	P2003-R2	P2003-L2	2 3/4" to 3"
Med	P2003-R3	P2003-L3	3" to 3 3/8"
Lrg	P2003-R4	P2003-L4	3 3/8" to 3 5/8"

3pp[®] Side Step[™] Splint

www.3pointproducts.com/3pp-side-step-splint/

Instructions: Measure length from web space to fingertip and width at proximal phalanx. Fits right or left hand. If between sizes, choose larger size.

	Right or Left	Web Space to Tip
Sm	P1201-2	1 1/2" - 2" (3.8-5cm)
Med	P1201-3	1 3/4" - 3 1/8" (4.4-8.3cm)
Lrg	P1201-4	3" - 3 1/2" (7.6-8.9cm)

Oval-8® Finger Splints

www.3pointproducts.com/oval-8-finger-splint--copy-1/

Please see the Oval-8 Size & Ordering Information sheet.

