

## Neb-u-Lite™ LX2 Nebulizer Compressor

- High-quality, low-noise aerosol compressor with user-friendly design.
- Piston-driven motor ensures high performance and long life.
- Complete with disposable nebulizer featuring 0.2 ml/min flow rate and 5 ml capacity medicine cup and cleaning/instruction manual.
- Designed to meet the requirements of HCPCS code: E0570.
- Five year limited warranty.

**JB0112-110** 1 ea



### Technical Specifications

Method Of Nebulization ..... Piston Driven  
 Power ..... 120VAC, 60Hz  
 Operating Temperature Range ..... 50°F TO 104°F (10°C TO 40°C)  
 Operating Humidity Range ..... 10% to 90% RH non-condensing  
 Storage Temperature Range ..... -5°F TO 122°F (-15°C TO 50°C)  
 Storage Humidity Range ..... 10% to 90% RH non-condensing  
 Weight ..... 4.2 lb (1.91 kg)

Width ..... 6.8"  
 Height ..... 4.2"  
 Length ..... 11.7"  
 Max Pressure (Compressor) ..... 30 psi  
 Flow Range (Compressor) ..... 0 - 14 liters/min  
 Particle Size (Nebulizer) ..... 0.5 to 5 microns  
 Capacity (Nebulizer) ..... 5 ml



**Graham-Field Health Products**

2935 Northeast Parkway  
 Atlanta, Georgia 30360

2007© Graham-Field Health Products

For Placing Orders:

**1-800-347-5678**

**Fax: 1-800-726-0601**

**www.grahamfield.com**

GF080001RevA08