



## Glass Dressing Jar 4" x 4"

4013

1/BX, 6 BX/CS

### Product Features

Provides exam room storage for products such as bandages, gauze, and cotton

Flint glass

Not autoclavable

Stainless steel lid

Each jar is individually boxed

### Attributes

<b>Brand Name</b>	Tech-Med
<b>Size</b>	4" x 4"
<b>Material</b>	Glass
<b>Sterile</b>	No
<b>Weight</b>	4.85 lbs

### Packaging

<b>Packaging</b>	1/BX, 6 BX/CS
<b>Inner Packaging</b>	1/BX
<b>Case Packaging</b>	6 BX/CS
<b>Total Case Quantity</b>	6